

[www.WARANTEC.com](http://www.WARANTEC.com)

**WARANTEC**

Inspiration from Pioneers

**IU/UT/KIT**  
**CATALOGUE** Ver.8.0



# WARANTEC


Inspiration from Pioneers

## THE FINE K-MED TECH

A DEEP COLLABORATION  
AMONG WORLD-CLASS  
PROFESSORS, CLINICIANS  
AND ENGINEERS







## NOTICE

1. Available products may vary by country due to registration process.
2. Images in this catalogue might be different from actual products.
3. Specifications are subject to change without prior notice.



# Contents

**HISTORY** 006

**INTRODUCTION** 007

**PRODUCT FEATURES** 008

## IU IMPLANT

**014**

IU  
Fixture



**023**

Cover  
Screw



**023**

Umbrella  
Cover Screw



## UT IMPLANT

**020**

UT  
Fixture



## PROSTHETICS

**024**

Healing  
Abutment



**027**

SB  
Abutment  
(Solid)



**028**

SB  
Abutment  
(Dual)



**031**

Milling  
Abutment



**032**

Angled  
Abutment



**034**

Temporary  
Abutment



**034**

Castable  
Abutment  
(CCM)



**035**

Impression  
Post  
Transfer  
Type



**035**

Impression  
Post  
Pick-up  
Type



**035**

Lab  
Analog



## InEx PROSTHETICS

IU/UT Flowchart

**038**

**040**

Single  
Multi-Unit  
Abutment



**041<sup>SB</sup>**

Abutment



**041**

Castable  
Abutment  
(CCM)



**041**

Temporary  
Cylinder



**042**

Impression  
Post Transfer  
Type



**042**

Impression  
Post Pick-Up  
Type



**042**

Lab Analog



**043**

Single-Unit  
Abutment



**044**

SB Abutment



**044**

Castable  
Abutment  
(CCM)



**044**

Temporary  
Cylinder



**045**

Impression  
Post Transfer  
Type



**045**

Impression  
Post Pick-Up  
Type



**045**

Lab Analog



**046**

Multi-Unit  
Abutment



**047**

Angled  
Multi-Unit  
Abutment



**048**

Milling  
Abutment



**048**

Multi-Unit  
SB  
Abutment



**049**

CCM Castable  
Abutment  
(CCM)



**049**

Temporary  
Cylinder



**050**

Impression  
Post Transfer  
Type



**050**

Impression  
Post Pick-Up  
Type



**050**

Lab Analog





**051**  
Adapter  
(Machine)



**051**  
Adapter  
(Torque Wrench)



## RETAINED

**053**  
Retained  
Abutment



## DIGITAL PROSTHETICS (Fixture level)

**055**  
Ti-Base  
Abutment



**057**  
Semicircle  
Ti-Base  
Abutment



**057**  
Custom  
Abutment



**058**  
Scan Body



**058**  
Digital Lab  
Analog



## DIGITAL PROSTHETICS (Abutment level)

**060**  
Ti-Base  
(Tapered  
Type)



**060**  
Ti-Base  
(Straight  
Type)



**061**  
Scan body



**061**  
Digital Lab  
Analog



IU Implant Surgical Procedure

IU Drill Sequence

UT Implant Surgical Procedure

UT Drill Sequence

IU Implant Surgical Procedure  
for Digital Dentistry

IU Drill Sequence for Digital  
Dentistry

IU Implant Surgical Procedure  
with Stopper Drill Kit

**062**

**067**

**071**

**075**

**079**

**082**

**087**

## SURGICAL KITS

**094**  
IU Full Kit



**104**  
IU  
Standard  
Kit



**106**  
IU  
Compact  
Kit



**108**  
IU  
Stopper Drill  
Kit

**118**  
IUT Kit



**122**  
UT  
Compact  
Kit



**126**  
WAGA Kit



## PROSTHETICS KITS

**137**  
WPK



**143**  
PM Kit



## OTHER KITS

**149**  
WISE Kit



**153**  
WISE Kit II



**160**  
WIRE Kit



**169**  
ESP Kit



**180**  
KAR Kit



## OTHER PRODUCTS

**192**  
WA Motor



**193**  
WA Bone



**94**  
Bone Dozer



**95**  
WA Mem



# HISTORY

## 1995~2000

---

G7 Government project for developing domestic implants

## 2001

---

Warantec was established based on the results of G7 Project

## 2002

---

IT system launched

## 2006

---

Obtained the patent for the ONEBODY system in the U.S  
(Patent No. US 6,981,873,B2)

## 2008

---

Obtained the ISO 13485, CE, and FDA certifications

## 2015

---

IU system launched

## 2017

---

YUHAN, No.1 Pharmaceutical company in Korea invested capital

## 2018

---

UT system launched

## 2020

---

KAR kit and ESP kit launched

## 2021

---

WISE kit, WIRE kit, IUT kit, and PM kit launched

## 2022

---

Obtained the patents for the KAR Kit and ESP Kit



# INTRODUCTION



Warantec, a leading dental implant company in Korea with a continued focus on R&D and clinical applications, was founded as a result of a 'next-generation artificial tooth development' research project supported by the Ministry of Health and Welfare as part of the Korean government's G7 Project in 1995 involving esteemed professors from the College of Dentistry of the Seoul National University, Yonsei University, and Catholic University.



Warantec's implants showcase outstanding and differentiated design, surface treatment (S.L.A.), and metal processing technology recognized by the implant academia since the beginning of its development. It spearheads the development of Korean implant design and surface treatment technologies.




Warantec introduced the first one-body ONEBODY IMPLANT SYSTEM for the first time in Korea through continued R&D investment and quality innovation. It has enhanced the success rate of implant procedures through its groundbreaking 7° connection IT IMPLANT SYSTEM. Furthermore, Warantec's IU IMPLANT SYSTEM, which has the 11° connection ensuring compatibility and convenience, has been acclaimed as an embodiment of breakthrough technologies.




Warantec has established successful clinical track record since 2001 and has been recognized for its value in the global market by attracting investments. These are clear testaments to Warantec's emerging presence in implantology at home and abroad.

## Inspiration from Pioneers

## Points to Implement When Designing our Implant



# Bone Collecting Drill





Designated bone collecting drill of Warantec



You can harvest enough autogenous bone for GBR while drilling.  
The best bone for the patient for free.




# GBR



The autogenous bone that is collected with the thanksgiving(bone collecting) drill makes the most stable ossification, and also can reduce the cost.



## Implant Stability




### Macro-design Concept

IU system has a combined design of straight and tapered, and UT system has micro-thread and progressive power thread design. Through these designs, each system can obtain the mechanically effective primary stability. After that, Warantec's certified SLA surface creates biologically stable osseointegration, and obtains quick secondary stability.


## Surface Treatment

Warantec


**S.L.A.**




### S.L.A. SURFACE




X500



X1,000



X5,000




X10,000

**S.L.A.** : Sand-blast **L**arge-grit **A**cid etching





# IU IMPLANT



## Characteristics of IU IMPLANTS

IU implant fixture has a design in which the straight-line thread and tapered body shape are properly combined.

By placing the maximum compression point 3mm below the cortical bone (subcortical area), the fixture can minimize the bone resorption of the cortical bone.

Also, by applying the square thread on the upper part of the fixture (cortical & subcortical area), the system can obtain excellent primary stability even on the soft bone, and decrease stress on the cortical bone.



### 01

Submerged type, internal hex and 11° tapered connection structure.

### 02

Disperse stress effectively by applying squared thread on the top so it can restrain the absorption of the alveolar bone.


















### 03

V-shaped thread at the bottom and optimized cutting-edge design for easier installation.

### 04

A body design to obtain good primary stability and proper stress distribution to prevent early bone resorption.

# IU Implant Specification

Color Coding	F3.3	F3.6	F4.0	F4.5	F5.0	F5.5	F6.0
							
	ORANGE	YELLOW	BLUE	PINK	PURPLE	GREEN	BROWN
							
UNIT : mm							
	3.1	3.4	3.8	4.3	4.8	5.3	5.8
Platform Diameter							
	3.3	3.6	4.1	4.6	5.1	5.6	6.1
Body Diameter							
	2.1	2.1	2.5	2.5	2.5	2.5	2.5
Semi Octa Deimension							




# IU Implant Fixture

• Cover screw is included in the package.

UNIT : mm


## F3.3 [MINI]

Art. No.	Length
FIU33085	8.5
FIU33100	10
FIU33115	11.5
FIU33130	13




## F3.6 [MINI]

Art. No.	Length
FIU36070	7.5
FIU36085	8.5
FIU36100	10
FIU36115	11.5
FIU36130	13




## F4.0

Art. No.	Length
FIU40060	7.5
FIU40070	7.5
FIU40085	8.5
FIU40100	10
FIU40115	11.5
FIU40130	13




## F4.5

Art. No.	Length
FIU45060	7.5
FIU45070	7.5
FIU45085	8.5
FIU45100	10
FIU45115	11.5
FIU45130	13




F5.0

Art. No.	Length
FIU50050	6.5
FIU50070	7.5
FIU50085	8.5
FIU50100	10
FIU50115	11.5
FIU50130	13



F5.5

Art. No.	Length
FIU55050	6.5
FIU55070	7.5
FIU55085	8.5
FIU50100	10
FIU55115	11.5
FIU55130	13



F6.0


Art. No.	Length
FIU60050	6.5
FIU60070	7.5
FIU60085	8.5
FIU60100	10







# UT IMPLANT



## Characteristics of **UT IMPLANTS**

UT implant system has a tapered body design of root form as a whole. The Progressive power thread that gets deepened from 2° to 4° of the lower apex makes an excellent primary stability and effectively disperses stress caused by the square shape. 400µm micro threads increase the secondary stability and disperse the stress by increasing the bone density and widening the surface area in the marginal bone.



### 01

A submerged type implant with an internal hex and 11° internal conical joint connection structure.

### 02

Disperse stress effectively by applying micro-thread so it can restrain the absorption of the alveolar bone.

### 03

A progressive power thread and square thread shape can disperse the stress and strengthen primary stability closer to the apex.

### 04


The root form shape and self-cutting edges increase the initial implantation stability and make the operation easier.

## UT Implant Specification

### Color Coding


F3.6	F4.3	F4.8	F5.3	F6.3
				
YELLOW	BLUE	PINK	PURPLE	BROWN
				

UNIT : mm



**Platform  
Diameter**

3.7	4.3	4.8	5.3	6.3
-----	-----	-----	-----	-----



**Bevel  
Height**

0.5	0.5	0.5	0.5	0.5
-----	-----	-----	-----	-----



**Hex  
Dimension**

2.1	2.5	2.5	2.5	2.5
-----	-----	-----	-----	-----


## UT Implant Fixture

• Cover screw is included in the package.

UNIT : mm


### F3.6 [MINI]

Art. No.	Length
FUT36070	7
FUT36085	8.5
FUT36100	10
FUT36115	11.5
FUT36130	13




### F4.3

Art. No.	Length
FUT43070	7
FUT43085	8.5
FUT43100	10
FUT43115	11.5
FUT43130	13




### F4.8

Art. No.	Length
FUT48070	7
FUT48085	8.5
FUT48100	10
FUT48115	11.5
FUT48130	13




### F5.3

Art. No.	Length
FUT53070	7
FUT53085	8.5
FUT53100	10
FUT53115	11.5
FUT53130	13



### F6.3



Art. No.	Length
FUT63070	7
FUT63085	8.5
FUT63100	10
FUT63115	11.5





## Common Platform for both IU and UT Implants

Both IU and UT Implants have 11-degree internal hexagonal connecting structure.



IU IMPLANT    UT IMPLANT



# Coding

Applicable Abutment and Shape

M

Mini

Hex

Non-Hex

Mini Hex

Mini Non-Hex

Mini - IU Implant : F3.3~F3.6

UT Implant : F3.6

S

Standard

Hex

Non-Hex

Standard Hex

Standard Non-Hex

Standard - IU Implant : F4.0~F6.0


UT Implant : F4.3~F6.3

22

## Cover Screw


- Assemble with 1.2 hex driver.
- Recommended tightening torque : under 10Ncm
- Package components : screw

Art. No.		Type
IUCSA		M1.6
IUCS		M2.0



## Umbrella Cover Screw

- Used in surgery accompanied by GBR and various bone augmentation procedures.
- Retaining the bone graft on the implant platform.
- Provides free direction according to the purpose of the practitioner.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : under 10Ncm
- Package components : Body + Screw




Art. No.	Diameter	G/H	Screw	Application
IUUCS5510N	5.5	1.0	UCS05	Standard
IUUCS5520N	5.5	2.0	UCS05	Standard
IUUCS5530N	5.5	3.0	UCS20	Standard
IUUCS6510N	6.5	1.0	UCS05	Standard
IUUCS6520N	6.5	2.0	UCS05	Standard
IUUCS6530N	6.5	3.0	UCS20	Standard




## Healing Abutment

- Assemble with 1.2 hex driver.
- Recommended tightening torque : under 10Ncm
- Package components : Abutment




### HA 4.0 [MINI] (D 3.95)

Art. No.	G/H
IUHA4010A	1.0
IUHA4020A	2.0
IUHA4030A	3.0
IUHA4040A	4.0
IUHA4050A	5.0
IUHA4060A	6.0



### HA 5.0 [MINI] (D 4.95)

Art. No.	G/H
IUHA5010A	1.0
IUHA5020A	2.0
IUHA5030A	3.0
IUHA5040A	4.0
IUHA5050A	5.0
IUHA5060A	6.0




**HA 4.0 (D 3.95)**

Art. No.	G/H
IUHA4010	1.0
IUHA4020	2.0
IUHA4030	3.0
IUHA4040	4.0
IUHA4050	5.0
IUHA4060	6.0

**HA 5.0 (D 4.95)**

Art. No.	G/H
IUHA5010	1.0
IUHA5020	2.0
IUHA5030	3.0
IUHA5040	4.0
IUHA5050	5.0
IUHA5060	6.0

**HA 6.0 (D 5.95)**


Art. No.	G/H
IUHA6010	1.0
IUHA6020	2.0
IUHA6030	3.0
IUHA6040	4.0
IUHA6050	5.0
IUHA6060	6.0

**HA 7.0 (D 6.95)**

Art. No.	G/H
IUHA7010	1.0
IUHA7020	2.0
IUHA7030	3.0
IUHA7040	4.0
IUHA7050	5.0
IUHA7060	6.0






## Solid / Dual



## SB Abutment (Solid)

- Abutment for the cement maintenance type prosthetics.
- Direct impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment




<b>SB 4.0</b> (D 3.95)	      <div> <div>M</div> <div> IUSB4015A IUSB4017A           </div> <div> IUSB4025A IUSB4027A           </div> <div> IUSB4035A IUSB4037A           </div> <div> IUSB4045A IUSB4047A           </div> <div> IUSB4055A IUSB4057A           </div> <div> IUSB4065A IUSB4067A           </div> </div>
<b>SB 4.0</b> (D 3.95)	      <div> <div>S</div> <div> IUSB4015 IUSB4017           </div> <div> IUSB4025 IUSB4027           </div> <div> IUSB4035 IUSB4037           </div> <div> IUSB4045 IUSB4047           </div> <div> IUSB4055 IUSB4057           </div> <div> IUSB4065 IUSB4067           </div> </div>
<b>SB 5.0</b> (D 4.95)	      <div> <div>S</div> <div> IUSB5015 IUSB5017           </div> <div> IUSB5025 IUSB5027           </div> <div> IUSB5035 IUSB5037           </div> <div> IUSB5045 IUSB5047           </div> <div> IUSB5055 IUSB5057           </div> <div> IUSB5065 IUSB5067           </div> </div>
<b>SB 6.0</b> (D 5.95)	      <div> <div>S</div> <div> IUSB6015 IUSB6017           </div> <div> IUSB6025 IUSB6027           </div> <div> IUSB6035 IUSB6037           </div> <div> IUSB6045 IUSB6047           </div> <div> IUSB6055 IUSB6057           </div> <div> IUSB6065 IUSB6067           </div> </div>



SB AButment (Dual)




































- Abutment for the cement & screw maintenance type prosthetics.
- Fixture level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package Components : Abutment + Screw




G/H A/H

<div>SB 4.0</div> <div>(D 3.95)</div> <div></div>	<div>IUSB4015HA IUSB4025HA IUSB4035HA IUSB4045HA IUSB4055HA IUSB4065HA IUSB4017HA IUSB4027HA IUSB4037HA IUSB4047HA IUSB4057HA IUSB4067HA</div>
<div>SB 4.0</div> <div>(D 3.95)</div> <div></div>	<div>IUSB4015NA IUSB4025NA IUSB4035NA IUSB4045NA IUSB4055NA IUSB4065NA IUSB4017NA IUSB4027NA IUSB4037NA IUSB4047NA IUSB4057NA IUSB4067NA</div>
<div>SB 4.0</div> <div>(D 3.95)</div> <div></div>	<div>IUSB4015H IUSB4025H IUSB4035H IUSB4045H IUSB4055H IUSB4065H IUSB4017H IUSB4027H IUSB4037H IUSB4047H IUSB4057H IUSB4067H</div>




<b>SB 4.0</b> (D 3.95)	       <div> IUSB4015N IUSB4025N IUSB4035N IUSB4045N IUSB4055N IUSB4065N  IUSB4017N IUSB4027N IUSB4037N IUSB4047N IUSB4057N IUSB4067N </div>
<b>SB 5.0</b> (D 4.95)	       <div> IUSB5015H IUSB5025H IUSB5035H IUSB5045H IUSB5055H IUSB5065H  IUSB5017H IUSB5027H IUSB5037H IUSB5047H IUSB5057H IUSB5067H </div>
<b>SB 5.0</b> (D 4.95)	       <div> IUSB5015N IUSB5025N IUSB5035N IUSB5045N IUSB5055N IUSB5065N  IUSB5017N IUSB5027N IUSB5037N IUSB5047N IUSB5057N IUSB5067N </div>
<b>SB 6.0</b> (D 5.95)	       <div> IUSB6015H IUSB6025H IUSB6035H IUSB6045H IUSB6055H IUSB6065H  IUSB6017H IUSB6027H IUSB6037H IUSB6047H IUSB6057H IUSB6067H </div>
<b>SB 6.0</b> (D 5.95)	       <div> IUSB6015N IUSB6025N IUSB6035N IUSB6045N IUSB6055N IUSB6065N  IUSB6017N IUSB6027N IUSB6037N IUSB6047N IUSB6057N IUSB6067N </div>

## Milling / Angled / Castable / Temporary



## Milling Abutment


- Abutment for the cement & screw maintenance type prosthetics.
- Easy margin formation.
- Fixture level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Screw



■ G/H ■ A/H

**MA 3.9**  
(D 3.9)

**MA 3.9**  
(D 3.9)



IUMA39LEA  
IUMA39LE

**MA 4.0**  
(D 3.95)

**MA 5.0**  
(D 4.95)

**MA 5.0**  
(D 4.95)

**MA 6.0**  
(D 5.95)



IUMA4017HA  
IUMA5017HA  
IUMA5017H  
IUMA6017H



IUMA4027HA  
IUMA5027HA  
IUMA5027H  
IUMA6027H



IUMA4037HA  
IUMA5037HA  
IUMA5037H  
IUMA6037H



IUMA4047HA  
IUMA5047HA  
IUMA5047H  
IUMA6047H




# Angled Abutment











- Abutment for the cement & screw maintenance type prosthetics.
- Compensates the fixture's implantation angle by 15° or 20°.
- Optional usage for various oral environments.
- Fixture level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Screw

G/H A/H

<p><b>AA 4.0_15°</b> (D 4.0)</p> 	<div>      </div> <div> <p>IUAA1526HA IUAA1536HA IUAA1546HA IUAA1556HA IUAA1566HA</p> <p>IUAA1528HA IUAA1538HA IUAA1548HA IUAA1558HA IUAA1568HA</p> </div>
<p><b>AA 4.0_20°</b> (D 4.0)</p> 	<div>      </div> <div> <p>IUAA2026HA IUAA2036HA IUAA2046HA IUAA2056HA IUAA2066HA</p> <p>IUAA2028HA IUAA2038HA IUAA2048HA IUAA2058HA IUAA2068HA</p> </div>








G/H A/H

<div>AA 5.0_15° (D 5.0)</div> <div>S</div>	<div></div> <div>IUAA1526H IUAA1528H IUAA1536H IUAA1538H IUAA1546H IUAA1548H IUAA1556H IUAA1558H IUAA1566H IUAA1568H</div>
<div>AA 5.0_20° (D 5.0)</div> <div>S</div>	<div></div> <div>IUAA2026H IUAA2028H IUAA2036H IUAA2038H IUAA2046H IUAA2048H IUAA2056H IUAA2058H IUAA2066H IUAA2068H</div>

## Temporary Abutment

- Abutment for the cement & screw maintenance type temporary prosthetics.
- Cut the cylinder when producing temporary prosthetics.
- Fixture level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw








Art. No.		Diameter	Type
IUTCHA		3.9	Hex
IUTCA		3.9	Non-Hex
IUTCH		3.9	Hex
IUTC		3.9	Non-Hex



## Castable Abutment (CCM)

- Abutment for the cement & screw maintenance prosthetics
- Castable for the nonprecious alloy, Use it for producing the customized prosthetics.
- Different colors per shape : Hex(Black), Non-hex(White)
- Fixture level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Packaged component : Abutment + Screw




Art. No.		Diameter	Type	Color
IUCA39CHA		3.9	Hex	Black
IUCA39CA		3.9	Non-Hex	White
IUCA39CH		3.9	Hex	Black
IUCA39C		3.9	Non-Hex	White



## Impression Post Transfer type (For Closed-Tray Technique)


- Cylinder shape makes it easier to transfer.
- Assemble with 1.2 hex driver.
- Package components : Impression Coping + Impression Screw



ART. NO.		Diameter	Length	Application	Included Screw
IUPTSA		4.0	11.6	IUHA4010A-IUHA4030A	IUSPTSA
IUPTLA		4.0	14.6	IUHA4040A-IUHA4060A	IUSPTLA
IUPTS40		4.0	11.6	IUHA4010-IUHA4030	IUSPTS40
IUPTL40		4.0	14.6	IUHA4040-IUHA4060	IUSPTL40
IUPTS		5.0	11.6	IUHA5010-IUHA5030	IUSPTS
IUPTL		5.0	14.6	IUHA5040-IUHA5060	IUSPTL

## Impression Post Pick-up type (For Open-Tray Technique)

- Rake shape for more precise and stable impression taking.
- Assemble by hand using 1.2 hex driver.
- Package components : Impression Coping + Impression Screw



ART. NO.		Length	Type	Included Screw
IUPPSA		12.2	Hex	IUSPPSA
IUPPLA		16.2	Hex	IUSPPLA
IUPPS		12.2	Hex	IUSPPS
IUPPL		16.2	Hex	IUSPPL


## Lab Analog

- Lab analog for the fixture level impression.
- Use it instead of a fixture when taking the impression.
- Optional usage according to the diameter of the fixture of F3.6 / F4.0 - F6.3.
- Package components : Lab analog

ART. NO.		Type
IUANM		Hex
IUANR		Hex



## InEx Prosthetics





# InEx Prosthetics

## Applicable Abutment and Coding



Hex



Non-Hex

**S4** stands for  
Single-unit Abutment  
D4.0



Hex

**S5** stands for  
Single-unit Abutment  
D5.0



Hex



Non-Hex

**M5** stands for  
Multi-unit Abutment  
D5.0











































Non-Hex

**AM5** stands for  
Angled Multi-unit Abutment  
D5.0


# IU / UT InEx Prosthetics Chart

		Single & Multi Unit	Impression	Analog		Tempor																																											
				Lab Analog	Digital Lab Analog																																												
 <b>SMUA 4.0</b>  Single & Multi Unit	Hex  Non Hex	 <table><tr><th>G/H</th><th>Art No.</th><th>G/H</th><th>Art No.</th></tr><tr><td>1.0</td><td>IUMU4010A</td><td>1.0</td><td>IUMU4010</td></tr><tr><td>2.0</td><td>IUMU4020A</td><td>2.0</td><td>IUMU4020</td></tr><tr><td>3.0</td><td>IUMU4030A</td><td>3.0</td><td>IUMU4030</td></tr><tr><td>4.0</td><td>IUMU4040A</td><td>4.0</td><td>IUMU4040</td></tr><tr><td>5.0</td><td>IUMU4050A</td><td>5.0</td><td>IUMU4050</td></tr><tr><td>6.0</td><td>IUMU4060A</td><td>6.0</td><td>IUMU4060</td></tr></table>	G/H	Art No.	G/H	Art No.	1.0	IUMU4010A	1.0	IUMU4010	2.0	IUMU4020A	2.0	IUMU4020	3.0	IUMU4030A	3.0	IUMU4030	4.0	IUMU4040A	4.0	IUMU4040	5.0	IUMU4050A	5.0	IUMU4050	6.0	IUMU4060A	6.0	IUMU4060	  Closed-Tray      Open-Tray <table><tr><th></th><th>Transfer Type</th><th>Pick-Up Type</th></tr><tr><td rowspan="3">Hex</td><td>ICOMUPT40SS</td><td>ICOMUPP40H</td></tr><tr><td>ICOMUPT40S</td><td></td></tr><tr><td>ICOMUPT40L</td><td></td></tr></table> Screw Included		Transfer Type	Pick-Up Type	Hex	ICOMUPT40SS	ICOMUPP40H	ICOMUPT40S		ICOMUPT40L		 ICOMAN40H  MUDAN40							
G/H	Art No.	G/H	Art No.																																														
1.0	IUMU4010A	1.0	IUMU4010																																														
2.0	IUMU4020A	2.0	IUMU4020																																														
3.0	IUMU4030A	3.0	IUMU4030																																														
4.0	IUMU4040A	4.0	IUMU4040																																														
5.0	IUMU4050A	5.0	IUMU4050																																														
6.0	IUMU4060A	6.0	IUMU4060																																														
	Transfer Type	Pick-Up Type																																															
Hex	ICOMUPT40SS	ICOMUPP40H																																															
	ICOMUPT40S																																																
	ICOMUPT40L																																																
 <b>SUA 5.0</b>  Single Unit	Hex	 <table><tr><th>G/H</th><th>Art No.</th></tr><tr><td>1.0</td><td>IUSU5010</td></tr><tr><td>2.0</td><td>IUSU5020</td></tr><tr><td>3.0</td><td>IUSU5030</td></tr><tr><td>4.0</td><td>IUSU5040</td></tr><tr><td>5.0</td><td>IUSU5050</td></tr><tr><td>6.0</td><td>IUSU5060</td></tr></table>	G/H	Art No.	1.0	IUSU5010	2.0	IUSU5020	3.0	IUSU5030	4.0	IUSU5040	5.0	IUSU5050	6.0	IUSU5060	  Closed-Tray      Open-Tray <table><tr><th></th><th>Transfer Type</th><th>Pick-Up Type</th></tr><tr><td rowspan="3">Hex</td><td>ICOMUPTSS</td><td>ICOMIPP43H</td></tr><tr><td>ICOMUPTS</td><td></td></tr><tr><td>ICOMUPTL</td><td></td></tr></table> Screw Included		Transfer Type	Pick-Up Type	Hex	ICOMUPTSS	ICOMIPP43H	ICOMUPTS		ICOMUPTL		 ICOSAN50H  SUDAN50																					
G/H	Art No.																																																
1.0	IUSU5010																																																
2.0	IUSU5020																																																
3.0	IUSU5030																																																
4.0	IUSU5040																																																
5.0	IUSU5050																																																
6.0	IUSU5060																																																
	Transfer Type	Pick-Up Type																																															
Hex	ICOMUPTSS	ICOMIPP43H																																															
	ICOMUPTS																																																
	ICOMUPTL																																																
 <b>MUA 5.0</b>  Multi Unit	Non Hex	 <table><tr><th>G/H</th><th>Art No.</th></tr><tr><td>1.0</td><td>IUMU5010</td></tr><tr><td>2.0</td><td>IUMU5020</td></tr><tr><td>3.0</td><td>IUMU5030</td></tr><tr><td>4.0</td><td>IUMU5040</td></tr><tr><td>5.0</td><td>IUMU5050</td></tr><tr><td>6.0</td><td>IUMU5060</td></tr></table>	G/H	Art No.	1.0	IUMU5010	2.0	IUMU5020	3.0	IUMU5030	4.0	IUMU5040	5.0	IUMU5050	6.0	IUMU5060	 Closed-Tray <table><tr><th></th><th>Transfer Type</th></tr><tr><td rowspan="3">Non-Hex</td><td>ICOMUPT43SS</td></tr><tr><td>ICOMUPT43S</td></tr><tr><td>ICOMUPT43L</td></tr></table> Screw Included		Transfer Type	Non-Hex	ICOMUPT43SS	ICOMUPT43S	ICOMUPT43L	 ICOMAN43  MUDAN50N																									
G/H		Art No.																																															
1.0	IUMU5010																																																
2.0	IUMU5020																																																
3.0	IUMU5030																																																
4.0	IUMU5040																																																
5.0	IUMU5050																																																
6.0	IUMU5060																																																
	Transfer Type																																																
Non-Hex	ICOMUPT43SS																																																
	ICOMUPT43S																																																
	ICOMUPT43L																																																
 <b>MUA 5.0</b>  Angled Multi-unit		 <table><tr><th colspan="2">Ø 5.0 (15°)</th><th colspan="2">Ø 5.0 (15°)</th></tr><tr><th>G/H</th><th>Art No.</th><th>G/H</th><th>Art No.</th></tr><tr><td>2.5</td><td>IUMUAA1625HA</td><td>2.5</td><td>IUMUAA1625H</td></tr><tr><td>3.5</td><td>IUMUAA1635HA</td><td>3.5</td><td>IUMUAA1635H</td></tr><tr><td>4.5</td><td>IUMUAA1645HA</td><td>4.5</td><td>IUMUAA1645H</td></tr></table> <table><tr><th colspan="2">Ø 5.0 (30°)</th><th colspan="2">Ø 5.0 (30°)</th></tr><tr><th>G/H</th><th>Art No.</th><th>G/H</th><th>Art No.</th></tr><tr><td>2.0</td><td>IUMUAA3030HA</td><td>2.0</td><td>IUMUAA3030H</td></tr><tr><td>4.0</td><td>IUMUAA3040HA</td><td>4.0</td><td>IUMUAA3040H</td></tr><tr><td>5.0</td><td>IUMUAA3050HA</td><td>5.0</td><td>IUMUAA3050H</td></tr></table>	Ø 5.0 (15°)		Ø 5.0 (15°)		G/H	Art No.	G/H	Art No.	2.5	IUMUAA1625HA	2.5	IUMUAA1625H	3.5	IUMUAA1635HA	3.5	IUMUAA1635H	4.5	IUMUAA1645HA	4.5	IUMUAA1645H	Ø 5.0 (30°)		Ø 5.0 (30°)		G/H	Art No.	G/H	Art No.	2.0	IUMUAA3030HA	2.0	IUMUAA3030H	4.0	IUMUAA3040HA	4.0	IUMUAA3040H	5.0	IUMUAA3050HA	5.0	IUMUAA3050H	 Open-Tray <table><tr><th></th><th>Pick-Up Type</th></tr><tr><td>Non-Hex</td><td>ICOMIPP43</td></tr></table> Screw Included		Pick-Up Type	Non-Hex	ICOMIPP43		
Ø 5.0 (15°)		Ø 5.0 (15°)																																															
G/H	Art No.	G/H	Art No.																																														
2.5	IUMUAA1625HA	2.5	IUMUAA1625H																																														
3.5	IUMUAA1635HA	3.5	IUMUAA1635H																																														
4.5	IUMUAA1645HA	4.5	IUMUAA1645H																																														
Ø 5.0 (30°)		Ø 5.0 (30°)																																															
G/H	Art No.	G/H	Art No.																																														
2.0	IUMUAA3030HA	2.0	IUMUAA3030H																																														
4.0	IUMUAA3040HA	4.0	IUMUAA3040H																																														
5.0	IUMUAA3050HA	5.0	IUMUAA3050H																																														
	Pick-Up Type																																																
Non-Hex	ICOMIPP43																																																

Y Cylinder	Milling & SB & CCM	Scan body	TI-base		Plastic cap	MUA grip	Adapter	Etc. Instrument
 Hex IOMUTC40H	<b>SB</b>  A/H      Hex 8        MUSB45H <b>CCM</b>  A/H      Hex IOMUCC40H	 MUSCT40S MUSCT40 MUSCT40L	Hex	Non-Hex	 ICOMUHC50	 ICOMUPC40	 HA-HHA  RA-HHA-S  RA-HHA-L  RA-HHA-LL	
			 MUTB45SHS MUTB45SH					
 Non-Hex IOMUTC40	<b>SB</b>  A/H      Non-Hex MUSB45N <b>CCM</b>  A/H      Non-Hex IOMUCC40			 MUTB45TNS MUTB45TN				
 Hex IOSUTC43H	<b>Milling</b>  Hex IOMUMA <b>SB</b>  Hex SUSB55H  IOSMUTC <b>CCM</b>  Hex IOSUCC50H	 MUSCT50S MUSCT50 MUSCT50L	 SUTB58SHS SUTB58SH					 WDE  RA-HS-S  RA-HS-L
 Non-Hex IOMUTC43	<b>Milling</b>  Non-Hex IOMUMAN <b>SB</b>  Non-Hex MUSB55N <b>CCM</b>  Non-Hex IOMUCC43	 MUSCT50NS MUSCT50N MUSCT50NL		 MUTB55TS MUTB55T	 ICOMUHC60	 ICOMUPC50	 HA-HH  RA-HH-S  RA-HH-L  RA-HH-LL	
				 MUTB55TL				
						 Grip : MUAAGR	 HD-TW-L	

## Single & Multi-Unit Abutment

- Abutment for the screw maintenance type prosthetics in a single case.
- Abutment level impression.
- Use the exclusive outer adapter for the single-unit abutment.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Grip + Healing Cap



G/H

**M SMUA4.0**  
(D 3.95)

**S SMUA4.0**  
(D 3.95)



IUMU4010A IUMU4020A IUMU4030A IUMU4040A IUMU4050A IUMU4060A  
IUMU4010 IUMU4020 IUMU4030 IUMU4040 IUMU4050 IUMU4060



Healing Cap : IC0MUHC50






Grip : IC0MUPC40

Healing Cap and Grip are included in the product package

## SB Abutment



- Abutment for the screw maintenance type prosthetics for the single-unit abutment.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw


ART. NO.		Diameter	Type
MUSB45H		4.5	Hex
MUSB45N		4.5	Non-Hex



## Castable Abutment (CCM)



- Abutment for the screw-retained type prosthetics for the single-unit abutment.
- Castable as nonprecious alloy, and it is used for producing the customized prosthetics.
- Added the abutment level impression & scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw


ART. NO.		Diameter	Type	Color
10MUCC40H		4.0	Hex	Black
10MUCC40		4.0	Non-Hex	White



## Temporary Cylinder

- Abutment for the screw-retained type temporary prosthetics for the single-unit abutment.
- Cut it to use when producing the temporary prosthetics.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw


ART. NO.		Diameter	Type
10MUTC40H		4.0	Hex
10MUTC40		4.0	Non-Hex








## Impression Post Transfer type (for Closed-Tray Technique)


- Transfer impression coping for single-unit abutment.
- Assemble with 1.2 hex driver.
- Package components : Impression Coping + Impression Screw




ART. NO.		Diameter	Length	Included Screw
ICOMUPT40SS		4.0	6.0	ISOMUPT40SS
ICOMUPT40S		4.0	8.5	ISOMUPT40S
ICOMUPT40L		4.0	11.5	ISOMUPT40L

## Impression Post Pick-Up Type (for Open-Tray Technique)


- Pick-up impression coping for single-unit abutment.
- Assemble with 1.2 hex driver.
- Package components : Impression Coping + Impression Screw \*
- Included Screw



ART. NO.		Diameter	Length	Type	Included Screw
ICOMUPP40H		4.9	7.7	Hex	ISOMIPPL

## Lab Analog


- Lab analog for single-unit Abutment level impression.
- Package components : Lab analog

ART. NO.		Diameter	Type
ICOMAN40H		4.0	Hex





## Single-Unit Abutment

- Abutment for the screw maintenance type prosthetics in a single case.
- Abutment level impression.
- Use the exclusive outer adapter for the single-unit abutment.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Grip + Healing Cap



 C/H

 **SUA 5.0**  
(D 4.95)



IUSU5010



IUSU5020



IUSU5030



IUSU5040



IUSU5050



IUSU5060



Healing Cap : IC0MUHC60





Grip : IC0MUPC50

Healing cap and grip are included in the product package

## SB Abutment


- Abutment for the screw maintenance type prosthetics for the single-unit abutment.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw


ART. NO.		Diameter	Type
SUSB55H		5.5	Hex



## Castable Abutment (CCM)

- Abutment for the screw-retained type prosthetics for the single-unit abutment.
- Castable as nonprecious alloy, and it is used for producing the customized prosthetics.
- Added the abutment level impression & scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw


ART. NO.		Diameter	Type	Color
IOSUCC50H		5.0	Hex	Black



## Temporary Cylinder


- Abutment for the screw-retained type temporary prosthetics for the single-unit abutment.
- Cut it to use when producing the temporary prosthetics.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw

ART. NO.		Diameter	Type
IO MUTC43H		5.0	Hex



## Impression Post Transfer type (for Closed-Tray Technique)


- Transfer impression coping for single-unit abutment.
- Assemble with 1.2 hex driver.
- Package components : Impression Coping + Impression Screw



ART. NO.		Diameter	Length	Included Screw
ICOMUPTSS	S5	5.0	6.0	ISOMUPTSS
ICOMUPTS	S5	5.0	8.5	ISOMUPTS
ICOMUPTL	S5	5.0	11.5	ISOMUPTL

## Impression Post Pick-Up Type (for Open-Tray Technique)

- Pick-up impression coping for single-unit abutment.
- Assemble with 1.2 hex driver.
- Package components : Impression Coping + Impression Screw \*
- Included Screw



ART. NO.		Diameter	Length	Type	Included Screw
ICOMIPP43H	S5	5.9	7.7	Hex	ISOMIPPL

## Lab Analog


- Lab analog for single-unit Abutment level impression.
- Package components : Lab analog

ART. NO.		Diameter	Type
IC0SAN50H	S5	5.0	Hex



## Multi-Unit Abutment

- Abutment for the screw maintenance type prosthetics in multiple cases.
- Abutment level impression.
- Assemble with the exclusive outer driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package Components : Abutment + Grip + Healing Cap



G/H

**S MUA5.0**  
(D 4.95)



IUMU5010



IUMU5020



IUMU5030



IUMU5040



IUMU5050



IUMU5060



Healing Cap : IC0MUHC60



Grip : IC0MUPC50

Healing Cap and Grip are included in the product package



## Angled Multi-Unit Abutment

- Abutment for the screw-retained type prosthetics of multiple cases.
- Same platform as the multi-unit abutment.
- Optional usage according to various oral environments. (16°, 30°)
- Abutment level impression & scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30 Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Grip + Screw + Healing Cap



Healing Cap : IC0MUHC60



Screw : IUMUAASA




IUMUAAS



Grip : MUAAGR

Screw and grip are included in the product package



G/H

**AMUA 5.0**  
(D 4.95)  
16°

**AMUA 5.0**  
(D 4.95)  
16°

IUMUAA1625HA  
IUMUAA1625HIUMUAA1635HA  
IUMUAA1635HIUMUAA1645HA  
IUMUAA1645H


**AMUA 5.0**  
(D 4.95)  
30°

**AMUA 5.0**  
(D 4.95)  
30°

IUMUAA3030HA  
IUMUAA3030HIUMUAA3040HA  
IUMUAA3040HIUMUAA3050HA  
IUMUAA3050H

## Milling Abutment


- Abutment for the screw maintenance type prosthetics for Multi-unit abutment.
- Easy margin formation.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Packaged component : Abutment + Screw



ART. NO.		Diameter	Type
IOMUMAN	M5 AM5	7.0	Non-Hex

## SB Abutment


- Abutment for the screw maintenance type prosthetics for Multi-unit abutment.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw



ART. NO.		Diameter	Type
MUSB55N	M5 AM5	5.5	Non-Hex

## Castable Abutment (CCM)


- Abutment for the screw-retained type prosthetics for the single-unit abutment.
- Castable as nonprecious alloy, and it is used for producing the customized prosthetics.
- Added the abutment level impression & scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw



ART. NO.		Diameter	Type	Color
IOMUCC43	M5	5.0	Non-Hex	White

## Temporary Cylinder


- Abutment for the screw-retained type temporary prosthetics for the single-unit abutment.
- Cut it to use when producing the temporary prosthetics.
- Abutment level impression.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw



ART. NO.		Diameter	Type
IOMUTC43	M5	5.0	Non-Hex

## Impression Post Transfer Type (for Closed-Tray Technique)


- Transfer impression coping for Multi-unit abutment.
- Assemble with 1.2 hex driver.
- Package Components : Impression Coping + Impression Screw




ART. NO.		Length	Type	Included Screw
ICOMUPT43SS	M5 AM5	6.0	Non-Hex	ISOMUPTSS
ICOMUPT43S	M5 AM5	5.9	Non-Hex	ISOMUPTS
ICOMUPT43L	M5 AM5	12	Non-Hex	ISOMUPTL

## Impression Post Pick-Up Type (for Open-Tray Technique)

- Pick-up impression coping for Multi-unit abutment.
- Assemble with 1.2 hex driver.
- Package Components : Impression Coping + Impression Screw




ART. NO.		Diameter	Type	Included Screw
ICOMIPP43	M5 AM5	5.9	Non-Hex	ISOMIPPL



## Lab Analog

- Lab analog for multi-unit abutment.
- Abutment level impression.
- Package Components : Lab Analog




ART. NO.		Diameter	Type
ICOMAN43	M5 AM5	5.0	Non-Hex




## InEx Instruments

### Adapter for Single-Unit / Multi-Unit Abutment (Machine)










- An exclusive outer adapter for single-unit / multi-unit abutments


ART. NO.		Diameter	Length
HA-HHA		4.0	22.3
HA-HH	 	5.0	22.3



### Adapter for Single-Unit / Multi-Unit Abutment (Torque Wrench)

- An exclusive outer hex adapter for single-unit / multi-unit abutments.
- Lengths of short, long, and extra long are available.

ART. NO.		Diameter	Length
RA-HHA-S		4.0	14.6
RA-HHA-L		4.0	19.1
RA-HHA-LL		4.0	24.0
RA-HH-S	 	5.0	14.6
RA-HH-L	 	5.0	19.1
RA-HH-LL	 	5.0	24.0






## Retained

If the housing and cap has the same specification with Zest's Locator, Rhein83's OT Equator Plus, and others, they are compatible with Warantec's product.



## Retained Abutment

- Abutment for the locator type overdenture.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30 Ncm after 10 minutes of initial placement.**
- Package components : Abutment



■ G/H

**M** RA4.3  
(D 4.3)

**S** RA4.3  
(D 4.3)



IULA4310AT  
IULA4310T

IULA4320AT  
IULA4320T


IULA4330AT  
IULA4330T

IULA4340AT  
IULA4340T

IULA4350AT  
IULA4350T

IULA4360AT  
IULA4360T

## Digital Prosthetics (Fixture Level)



### Applicable software


3shape  
exocad

Downloadable on  
warantec.com



## Ti-Base Abutment


- Abutment for the cement & screw maintenance type prosthetics.
- Used for producing the custom abutment with the CAD / CAM devices.
- Fixture level scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package Components : Abutment + Screw



G/H

**M** TB4.0  
(D 4.0)

**S** TB4.0  
(D 4.0)



IUTB4015HA IUTB4025HA IUTB4035HA IUTB4045HA IUTB4055HA IUTB4065HA  
IUTB4015H IUTB4025H IUTB4035H IUTB4045H IUTB4055H IUTB4065H


**M** TB4.0  
(D 4.0)


**S** TB4.0  
(D 4.0)




IUTB4015NA IUTB4025NA IUTB4035NA IUTB4045NA IUTB4055NA IUTB4065NA  
IUTB4015N IUTB4025N IUTB4035N IUTB4045N IUTB4055N IUTB4065N

G/H

 **TB5.0**  
(D 5.0)  
Hex




IUTB5015H   IUTB5025H   IUTB5035H   IUTB5045H   IUTB5055H   IUTB5065H

 **TB5.0**  
(D 5.0)  
Non-Hex




IUTB5015N   IUTB5025N   IUTB5035N   IUTB5045N   IUTB5055N   IUTB5065N

 **TB5.8**  
(D 5.8)  
Hex



IUTB5815H   IUTB5825H   IUTB5835H   IUTB5845H   IUTB5855H   IUTB5865H


 **TB5.8**  
(D 5.8)  
Non-Hex



IUTB5815N   IUTB5825N   IUTB5835N   IUTB5845N   IUTB5855N   IUTB5865N

## Semicircle Ti-Base Abutment


- Abutment for the cement & screw maintenance type prosthetics.
- Fixture level scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Screw



Art. No.		G/H
IUTB4510HA		1
IUTB4525HA		2.5
IUTB4510H		1
IUTB4525H		2.5

## Custom Abutment

- Abutment for the cement & screw maintenance type prosthetics.
- Used for producing the custom abutments with the CAD / CAM devices.
- Extra abutment screws are included for lab.
- Fixture level scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 30Ncm
- **Re-tighten with a torque of 30Ncm after 10 minutes of initial placement.**
- Package components : Abutment + Screw x 2ea








Art. No.		Art. No.		Diameter	Machine Type
IUCMAA		IUCMA		10	Arum
IUCMA14A		IUCMA14		14	Arum
IUCMBA		IUCMB		10	Manix
IUCMCA		IUCMC		12	Medentika
IUCMDA		IUCMD		14	Denkai







## Scan Body


- Scan body for the digital abutment / digital crown.
  - Application for the chair side / lab side.
  - Assemble with 1.2 hex driver using hand force.
  - Fixture level scan.
  - Package components : Scan body + Screw
- Since there is a horizontal groove in the middle of the scan body, conventional impression-taking method is available.



Art. No.		Diameter	Length	Application
IUSCTHA		4.3	8.5	Short
IUSCTHLA		4.0	13.7	Long
IUSCTH		4.3	8.5	Short
IUSCTHL		4.0	13.7	Long



## Scan Body Screw


Art. No.		Length	Application
IUSCSA		14.4	Short
IUSCLA		19.6	Long
IUSCS		14.4	Short
IUSCL		19.6	Long




## Digital Lab Analog

- Digital lab analog for printed model.
- Package components : Analog + Screw

Art. No.		Diameter	Length
IUDANA		4.5	9
IUDAN		4.5	9




## InEx Digital Prosthetics (Abutment Level)



## Ti-Base (Tapered Type)

- Abutment for the cement & screw maintenance type prosthetics.
- Used for the digital abutments with CAD / CAM devices.
- Abutment level scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw




ART. NO.		Diameter(mm)	A/H
MUTB45TNS	S4	4.5	4.0
MUTB45TN	S4	4.5	5.5
MUTB55TS	M5 AM5	5.5	4.2
MUTB55T	M5 AM5	5.5	6.0
MUTB55TL	M5 AM5	5.5	11.0



## Ti-Base (Straight Type)

- Abutment for the cement & screw maintenance type prosthetics.
- Used for the digital abutment with the CAD / CAM devices.
- Abutment level scan.
- Assemble with 1.2 hex driver.
- Recommended tightening torque : 20Ncm
- Package components : Abutment + Screw




ART. NO.		Diameter(mm)	A/H
MUTB45SHS	S4	4.5	4.0
MUTB45SH	S4	4.5	5.8
SUTB58SHS	S5	5.8	4.2
SUTB58SH	S5	5.8	6.3

## Scan Body

- Scan body for digital prosthetics.
- Joint usage for Single-Unit Abutment / Multi-Unit Abutment.
- Application for the chair side / lab side.
- Abutment level scan.
- Assemble with 1.2 hex driver.
- Package components : Scan Body + Screw

Art. No.		Diameter	Length
MUSCT40S	S4	4.0	5.0
MUSCT40	S4	4.0	7.0
MUSCT40L	S4	4.0	11
MUSCT50S	S5	5.0	5.0
MUSCT50	S5	5.0	7.0
MUSCT50L	S5	5.0	11
MUSCT50NS	M5 AM5	5.0	5.0
MUSCT50N	M5 AM5	5.0	7.0
MUSCT50NL	M5 AM5	5.0	11



## Digital Lab Analog

- Digital lab analog for abutment level scan.
- Package components : Analog + Screw



Art. No.		Diameter
MUDAN40	S4	4.0
MUDAN50N	M5 AM5	5.0
SUDAN50	S5	5.0

# IU IMPLANT SURGICAL PROCEDURE


It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

# IU Implant Surgical Procedure (In Case of Fixture 4.5x10mm)

Please follow this procedure for a successful and safe operation.


## 1. Point Drill

- Recommended RPM : 800-1,000 RPM
- Mark the implantation place for the fixture on the cortical bone.  
(In Case of Fixture 4.5X10mm)




## 2. Lindemann Drill

- Drill it to the laser-marked line equal to the fixture length.
- Use a longer drill or extension if adjacent tooth blocks the drill path.
- Use the drill at 1,000 RPM, and increase the amount of irrigation to reduce bone overheating as high RPM causes more heat.
- Angle correction is available while drilling.



## 3. Guide Pin


- Used to identify the direction and placement.
- Suggestion : insert dental floss into the hole in the middle to stably fix the pin.






#### 4. Pilot Drill

- Suggestion : 800-1,000 RPM
- Drill for extending the hole for the next step.
- Maintain the drilling axis when using a drill.
- Check the laser mark on the cortical bone and use it intuitively.
- Can be used on the hard bone selectively.




#### 5. Bone Collecting Drill F3.6 X 10mm

- Suggestion : 800-1,000 RPM
- Drill it to the laser marked line equal to the fixture length as an intermediate drilling process.
- Select the length of short or long specifications according to the intraoral access environment.
- It is recommended to use a stopper for more bone harvest from irrigation.




#### 6. Bone Collecting Drill F4.0 X 10mm

- Suggestion : 800-1,000 RPM
- Drill it to the laser marked line equal to the fixture length as an intermediate drilling process.
- Select the length of short or long specifications according to the intraoral access environment.
- It is recommended to use a stopper for more bone harvest from irrigation.




## 7. Bone Collecting Drill F4.5 X 10mm

- Suggestion : 800-1,000 RPM
- Select the length of short or long specifications according to the intraoral access environment.
- Drill it to the laser marked line equal to the fixture length as an intermediate drilling process.
- Select the length of short or long specifications according to the intraoral access environment.
- It is recommended to use a stopper for more bone harvest from irrigation.



## 8. Final Drill

- Suggestion : 800-1,000 RPM
- Drill for cortical bone's hole expansion after using the bone collecting drill.
- Use it on hard bones only.
- Use designated fixtures by diameter.




## 9. Move the fixture from the ampule into the oral cavity

- When the user opens the lid of the ampule, the implant gets exposed and is ready to be grabbed.
- As it is a no-mount type fixture, use the fixture driver to grab the fixture firmly by feeling the hexes are met, and lift it from the ampule.
- Face the fixture upwards to prevent dropping it when moving the fixture in to the oral cavity.




## 10. Placement of the F4.5 X 10mm fixture

- Can select the machine type or torque wrench type fixture driver to implant the fixture.
- Recommended implantation torque : 25 - 40 Ncm



## 11. Tightening Cover Screw

- Open the bottom lid of the ampule and grab the cover screw with 1.2 hex driver.
- Face the cover screw upwards to prevent dropping it when moving the cover screw into the oral cavity.
- Tighten the cover screw with a torque of 10 Ncm or less.




## Surgical Procedure for IU IMPLANT

After using Lindemann Drill, the operator can use Pilot Drill (PD-30) selectively.


It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

### IU F3.3 Surgical Sequence



Bone Quality	Point	Lindemann	Straight (F3.3)	Final (F3.3)	Fixture Driver
Hard	▶	▶	▶	▶	Implant Placement
Normal		▶	▶		
Soft		▶	▷		


### IU F3.6 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Final (F3.6)	Fixture Driver
Hard	▶	▶	▶	▶	Implant Placement
Normal		▶	▶		
Soft		▶	▷		


It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

IU F4.0 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Final (F4.0)	Fixture Driver
Hard	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶		
Soft		▶	▶			


IU F4.5 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Collecting (F4.5)	Final (F4.5)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶		
Soft		▶	▶	▶			


It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

## IU F5.0 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Collecting (F4.5)	Collecting (F5.0)	Final (F5.0)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶	▶		
Soft		▶	▶	▶	▶			

## IU F5.5 Surgical Sequence




Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Collecting (F4.5)	Collecting (F5.0)	Collecting (F5.5)	Final (F5.5)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶	▶	▶		
Soft		▶	▶	▶	▶	▶			



It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

IU **F6.0** Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Collecting (F4.5)	Collecting (F5.0)	Collecting (F5.5)	Collecting (F6.0)	Final (F6.0)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶	▶	▶	▶		
Soft		▶	▶	▶	▶	▶	▶			



# **UT IMPLANT SURGICAL PROCEDURE**

It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

# UT Implant Surgical Procedure (In Case of Fixture 4.3x10mm)

Please follow this procedure for a successful and safe operation.


## 1. Point Drill

- Suggestion : 800-1,000 RPM
- Mark the implantation place for the fixture on the cortical bone.




## 2. Lindemann Drill

- Drill it to the laser-marked line equal to the fixture length.
- Use a longer drill or extension if adjacent tooth blocks the drill path.
- Use the drill at 1,000 RPM, and increase the amount of irrigation to reduce bone overheating high RPM causes more heat.
- Angle correction is available while drilling.




## 3. Guide Pin

- Use it in order to identify the direction and placement.
- Suggestion : insert a floss into the hole in the middle to fix the pin stably.




#### 4. Bone Collecting Drill F3.6 X 10mm

- Suggestion : 800-1,000 RPM
- Drill it to the laser marked line equal to the fixture length as an intermediate drilling process.
- Select the length of short or long specifications according to the intraoral access environment.
- It is recommended to use a stopper for more bone harvest from irrigation.




#### 5. Final Drill F4.0

- Suggestion : 800-1,000 RPM
- Drill for cortical bone expansion after using the bone collecting drill.
- Use it on hard bones only.
- Use designated fixtures by diameter.



#### 6. Tap Drill F4.3 x 7mm

- Tap for tapered fixtures.
- It is used to create a micro threaded channel in hard bone.
- Use machine or torque wrench assembled with shank adapter to tap it.




## 8. Move the fixture from the ampule into the oral cavity

- When the user opens the lid of the ampule, the implant gets exposed and is ready to be grabbed.
- As it is a no-mount type fixture, use the fixture driver (magic grip) to grab the fixture and lift it from the ampule.
- Face the fixture upwards to prevent dropping it when moving the fixture in to the oral cavity.




## 9. Placement of the F4.3 X 10mm fixture

- Can select the machine type or torque wrench type fixture driver to implant the fixture.
- Recommended implantation torque : 25 - 40 Ncm



## 10. Tightening Cover Screw


- Open the bottom lid of the ampule and grab the cover screw with 1.2 hex driver.
- Face the cover screw upwards to prevent dropping it when moving the cover screw into the oral cavity.
- Tighten the cover screw with a torque of 10Ncm or less.



## Surgical Procedure for UT IMPLANT


It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

### UT F3.6 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Final (F3.6)	Tap (F3.6)	Fixture Driver
Hard	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶		
Soft		▶	▶			


### UT F4.3 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Final (F4.0)	Tap (F4.3)	Fixture Driver
Hard	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶		
Soft		▶	▶			


It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

UT F4.8 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Final (F4.5)	Tap (F4.8)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶		
Soft		▶	▶	▶			


UT F5.3 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Collecting (F4.5)	Final (F5.0)	Tap (F5.3)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶	▶		
Soft		▶	▶	▶	▶			

It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

### UT F6.3 Surgical Sequence



Bone Quality	Point	Lindemann	Collecting (F3.6)	Collecting (F4.0)	Collecting (F4.5)	Collecting (F5.0)	Collecting (F5.5)	Final (F6.0)	Tap (F6.3)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal		▶	▶	▶	▶	▶	▶	▶		
Soft		▶	▶	▶	▶	▶	▶			






# **IU IMPLANT SURGICAL PROCEDURE FOR DIGITAL DENTISTRY**

## IU Surgical Procedure of IU for Digital Dentistry


- 2 types of metal sleeves are provided according to the diameter of the guide hole.
  - Regular Sleeve : F3.6/F4.0/F4.5/F5.0
  - Wide Sleeve : F5.5/F6.0
- Select the length of drills according to the offset.

		Fixture Length (mm)				
		7.5	8.5	10	11.5	13
Offset (mm)		↓	↓	↓	↓	↓
10	7.5	Used Drill Length : Use a drill with the same length as the fixture.				
	7.5	8.5	10	11.5	13	
11.5	8.5	Used Drill Length : Use a drill 1.5mm longer than the fixture length.				
	8.5	10	11.5	13	14.5	
13	10	Used Drill Length : Use a drill 3mm longer than the fixture length.				
	10	11.5	13	14.5	14.5 (Drill it up to 16)	

- The L7.5mm fixture has offsets of 10/11.5/13mm, distinguishing it from other lengths of fixtures. It shares the same surgical procedure as the F8.5 fixture.




- 14.5mm drill does not have a drill stopper, and can drill up to 16mm by using the 1.5mm marked line.




- The drills that are included in the WAGA kit has a standard length of Offset 9mm. Therefore, use the drill selectively in accordance with the guide and offset below.
- Warantec Guide can use sleeve and sleeveless selectively.


### Offset 9mm



### Offset 10.5mm




### Offset 12mm




## IU F3.6X10mm Surgical Sequence

### Offset 9mm




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 10mm	Final Drill F3.6	Fixture Driver
Hard	►	►	►	►	►	Implant Placement
Normal	►	►	►	►		
Soft	►	►	►	►		

### Offset 10.5mm




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 11.5mm	Final Drill F3.6	Fixture Driver
Hard	►	►	►	►	►	Implant Placement
Normal	►	►	►	►		
Soft	►	►	►	►		

### Offset 12mm




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 13mm	Final Drill F3.6	Fixture Driver
Hard	►	►	►	►	►	Implant Placement
Normal	►	►	►	►		
Soft	►	►	►	►		


## IU F4.0X10mm Surgical Sequence



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 10mm	F4.0 10mm	Final Drill F4.0	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶			




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 11.5mm	F4.0 11.5mm	Final Drill F4.0	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶			




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 13mm	F4.0 13mm	Final Drill F4.0	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶			


## IU F4.5X10mm Surgical Sequence



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 10mm	F4.0 10mm	F4.5 10mm	Final Drill F4.5	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	
Soft	▶	▶	▶	▶	▶	▶	▶	




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 11.5mm	F4.0 11.5mm	F4.5 11.5mm	Final Drill F4.5	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	
Soft	▶	▶	▶	▶	▶	▶	▶	




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 13mm	F4.0 13mm	F4.5 13mm	Final Drill F4.5	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	
Soft	▶	▶	▶	▶	▶	▶	▶	




## IU F5.0X10mm Surgical Sequence



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 10mm	F4.0 10mm	F4.5 10mm	F5.0 10mm	Final Drill F5.0	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	▶	
Soft	▶	▶	▶	▶	▶	▶	▶	▶	



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 11.5mm	F4.0 11.5mm	F4.5 11.5mm	F5.0 11.5mm	Final Drill F5.0	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	▶	
Soft	▶	▶	▶	▶	▶	▶	▶	▶	



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	F3.6 13mm	F4.0 13mm	F4.5 13mm	F5.0 13mm	Final Drill F5.0	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	▶	
Soft	▶	▶	▶	▶	▶	▶	▶	▶	






# **IU IMPLANT SURGICAL PROCEDURE WITH STOPPER DRILL KIT**

Surgical Sequence for IU IMPLANT

IU Stopper Drill




IU F3.3 Surgical Sequence




Bone Quality	NLD-20	NSD-20 (ø2.7)	NPD-33 (F3.6)	Fixture Driver
Hard	▶	▶	▶	Implant Placement
Normal	▶	▶		
Soft	▶			

### IU F3.6 Surgical Sequence




Bone Quality	Point Drill	Initial Drill (ø2.2)	Final Drill (F3.6)	Cortical Drill (F3.6)	Fixture Driver
Hard	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶		
Soft	▶	▶			

### IU F4.0 Surgical Sequence




Bone Quality	Point Drill	Initial Drill (ø2.2)	Final Drill (F3.6)	Final Drill (F4.0)	Cortical Drill (F4.0)	Fixture Driver
Hard	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶		
Soft	▶	▶	▶			

## IU F4.5 Surgical Sequence




Bone Quality	Point Drill	Initial Drill (ø2.2)	Final Drill (F3.6)	Final Drill (F4.0)	Final Drill (F4.5)	Cortical Drill (F4.5)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶			

## IU F5.0 Surgical Sequence




Bone Quality	Point Drill	Initial Drill (ø2.2)	Final Drill (F3.6)	Final Drill (F4.0)	Final Drill (F4.5)	Final Drill (F5.0)	Cortical Drill (F5.0)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶	▶			

## IU F5.5 Surgical Sequence



Bone Quality	Point Drill	Initial Drill (ø2.2)	Final Drill (F3.6)	Final Drill (F4.0)	Final Drill (F4.5)	Final Drill (F5.0)	Final Drill (F5.5)	Cortical Drill (F5.5)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶	▶	▶			

## IU F6.0 Surgical Sequence



Bone Quality	Point Drill	Initial Drill (ø2.2)	Final Drill (F3.6)	Final Drill (F4.0)	Final Drill (F4.5)	Final Drill (F5.0)	Final Drill (F5.5)	Collecting (F6.0)	Final (F6.0)	Fixture Driver
Hard	▶	▶	▶	▶	▶	▶	▶	▶	▶	Implant Placement
Normal	▶	▶	▶	▶	▶	▶	▶	▶		
Soft	▶	▶	▶	▶	▶	▶	▶			





# Kit

---

## **SURGICAL KIT**

- 1. IU FULL KIT**
- 2. IU STANDARD KIT**
- 3. IU COMPACT KIT**
- 4. IU STOPPER DRILL KIT**
- 5. IUT KIT**
- 6. UT COMPACT KIT**
- 7. WAGA KIT**

---

## **PROSTHETICS KIT**

- 8. WPK**  
(WARANTEC PROSTHETICS KIT)
- 9. PM KIT**  
(PROSTHETICS MULTI KIT)

---

## **OTHERS**

- 10. WISE KIT**  
(WARANTEC IMPLANT'S SINUS ELEVATION KIT)
  - 11. WISE KIT II**  
(WARANTEC IMPLANT'S SINUS ELEVATION KIT II)
  - 12. WIRE KIT**  
(WARANTEC IMPLANT'S RESCUE FOR EMERGENCY KIT)
  - 13. ESP KIT**  
(EASY SCREWMENTATION POSITIONER KIT)
  - 14. KAR KIT**  
(KINETIC ABUTMENT REMOVER KIT)
-



# IU Full Kit


## Full Kit for IU Implant




It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

# Full kit for IU implant(F3.6-F6.0)

Art No. : IU KIT



1 Point Drill Lindemann Drill Pilot Drill



2 LD & BC Drill Stopper



3 Bone Collecting Drill



4 Final Drill




5 G/H Checker & Path Finder



6 Guide Pin



7 Fixture Driver



8 Final Drill Stopper




9 Extension




10 1.2 Hex



11 Torque Wrench




(Bottom of the Case)  
Depth Gauge



Point Drill

- Initial drill for marking the implantation point easily on the cortical bone without slipping when making an implantation hole.


Art No.	Diameter(mm)	Edge Length(mm)
PD	1.5	9.5



Lindemann Drill

- Make a drill hole to determine the implantation depth and axis.
- Drilling as deep as the length of the fixture is required.
- Short and long drill are provided for the convenience of the user.
- Angle correction is available while drilling.


Art No.	Length(mm)	Type
LD-20-S	30	Short
LD-20-L	33	Long



Pilot Drill

- Maintain the path of the previous drill when using the next drill.
- There is a guide at the bottom of the drill to widen the diameter of the upper part of the hole with a small diameter, and the top of the drill has a straight blade wider than the guide.


Art No.
PD-30



## Bone Collecting Drill

- Sufficient amount of autogenous bone can be harvested simultaneously when drilling.
- Designed to minimize the heat while drilling.
- It is easy to adjust the drilling depth by using the drill stopper.
- Short and long drills are provided for the convenience of the user.

Fixture	F3.6	F4.0	F4.5	F5.0	F5.5	F6.0
Short	SD-30-S	SD-35-S	SD-40-S	SD-45-S	SD-50-S	SD-55-S
Long	SD-30-L	SD-35-L	SD-40-L	SD-45-L	SD-50-L	SD-55-L
Diameter(mm)	3.0	3.6	4.1	4.6	5.2	5.7



## Final Drill


- Drill for expanding the cortical bone after using the bone collecting drill.
- Use it for making the final holes on hard bones.
- Use in consideration of the diameter of the fixture.



Fixture	F3.6	F4.0	F4.5	F5.0	F5.5	F6.0
Art No.	PD-35	PD-40	PD-45	PD-50	PD-55	PD-60
Diameter(mm)	3.6	4.1	4.6	5.1	5.6	6.1

Guide Pin

- Used to identify the direction and placement of the bone preparation after using the Lindemann Drill.



Art No.
GPS

Drill Extension

- An instrument to extend the length of the drill and the driver of the machine.
- Incomplete assembly may cause bending or fracture of the extension.
- The length can be extended by 15.4mm when the drill extension is assembled.




Art No.	Extendable Length(mm)
WDE	15.4

Gingiva Height Checker & Path Finder

- Check the height of the gingiva with the indicated color line.
- In the multiple case, assemble with the fixture that is placed firstly and confirm the precise implantation path.

Art No.	Application
IC0GHCA	MINI
IC0GHC	STANDARD



## Bone Collecting Drill Stopper

- The length of the drill stopper indicates the actual length remaining when the drill stopper is assembled with the bone collecting drill.
- Color coding is applied to each stopper for the convenience of users to easily identify the length and relocate it back into the kit.



Art No.							Application
DS-5	DS-6	DS-7	DS-8.5	DS-10	DS-11.5	DS-13	F3.3-F5.0
DSW-5	-	DSW-7	DSW-8.5	DSW-10	DSW-11.5	DSW-13	F5.5-F6.0

Optional

- This stopper makes actual length stop which excluding the Y-dimension length.
- It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

## Final Drill Stopper

- The length of the drill stopper indicates the depth of bone to be removed when the drill stopper is assembled with the final drill.




Art No.	Application
PS-3.5	F3.3-F5.0
PSW-3.5	F5.5-F6.0

## Fixture Driver Machine Type


- Used when placing the fixture with a machine.
- Connections are divided into F3.3(Mini)/F4.0-6.0(Standard).
- Short and long drivers are provided for the convenience of the user.

Art No.	Length(mm)	Application	
HD-IU-MS	25	Mini	* Optional
HD-IU-ML	29		
HD-IU-RS	25	Standard	* Optional
HD-IU-RL	29		



## Fixture Driver Torque Wrench Type

- Used when placing the fixture with a torque wrench.
- Must check that it is completely assembled before use.  
Incomplete assembly may cause fracture due to heavy stress.
- Be careful, as excessive torque may cause damage to the internal hex.
- Short and long drivers are provided for the convenience of the user.




Art No.	Length(mm)	Application	
RD-IU-MS	13.2	Mini	* Optional
RD-IU-ML	18.8		
RD-IU-RS	13.2	Standard	* Optional
RD-IU-RL	18.8		

Depth Gauge

- Instrument for measuring the drilling depth. (7-13mm)


Art No.
DP



Torque Wrench

- Assemble it to the adapter or driver to check torque during implantation or screw tightening.
- Apply torque by pulling the bar to the indicated torque value line the user wants to apply.
- Rotate the pivot hand of the torque wrench to change the torque direction.

Art No.
TW






## 1.2 Hex Driver Machine Type

\* Optional


- Driver for Machine.
- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 45Ncm. Exceeding the torque range may damage the tip.
- Various lengths are available for the convenience of the user.

Art.No	Length(mm)
HD-LT-S	20
HD-LT-L	24
HD-LT-LL	28



## 1.2 Hex Driver Torque Wrench type

- Can be used for both torque wrench & hand driver.
- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 45 Ncm. Exceeding the torque range may damage the tip.
- Various lengths are available for the convenience of the user.



Art.No	Length(mm)	
HD-TW-S	8mm	* Optional
HD-TW-L	13mm	
HD-TW-LL	18mm	* Optional


## ADDITIONAL PRODUCTS

## Narrow Lindemann Drill

\* Optional

- Lindemann drill designated for IU F3.3.
- Make a drill hole to determine the implantation depth and axis.
- Drilling as deep as the length of the fixture is required.
- Angle correction is available while drilling.

Art.No	Application
NLD-20	IU F3.3




## Narrow Straight Drill

\* Optional

- Straight drill designated for IU F3.3.
- Drilling as deep as the length of the fixture is required.
- The color coding on the upper part of drill is indicating the diameter and the main fixture that is being used.

Art.No	Application
NSD-27	IU F3.3




## Narrow Final Drill

\* Optional

- Final drill designated for IU F3.3.
- Drill for expanding the cortical bone.
- Use it after making the final holes on hard bones.

Art.No	Application
NPD-33	IU F3.3



# IU Standard Kit


Standard Kit for IU Implant



It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

# Standard Kit for IU Implant(F3.6-F5.0)

Art No. : IU KIT-S



1 Lindemann Drill



LD-20-S

2 Pilot Drill




PD-30

3 Bone Collecting Drill



F3.6 SD-30-S F4.0 SD-35-S F4.5 SD-40-S F5.0 SD-45-S

4 Final Drill



PD-35 PD-40 PD-45 PD-50

5 LD & BC Drill Stopper




DS-7 DS-8.5 DS-10 DS-11.5 DS-13

6 Pilot/Final Drill Stopper



PS-3.5

7 Fixture Driver




HD-IU-ML RD-IU-ML HD-IU-RL RD-IU-RL RD-IU-RL

8 Guide Pin



GPS

9 G/H Checker & Path Finder



IC0GHCA IC0GHC

10 1.2 Hex




HD-TW-L

11 Extension



WDE



(Bottom of the Case)  
Torque Wrench

Depth Gauge

TW

DP

# IU Compact Kit


## Compact Kit for IU Implant



It is always recommended to use 1 step shorter stopper for drilling than the length of fixture to be installed, to make sure the fixture will be submerged 0.5~1.0mm deeper than the bone level.

# Compact Kit for IU Implant(F4.0-F5.0)

Art No. : IU BASIC KIT



## 1 Lindemann Drill



LD-20-S

## 2 Pilot Drill



PD-30

## 3 Bone Collecting Drill



SD-35-S



SD-40-S



SD-45-S

## 4 Final Drill



PD-40



PD-45



PD-50

## 5 Guide Pin



GPS

## 6 G/H Checker & Path Finder



ICOGHC

## 7 LD & BC Drill Stopper



DS-7



DS-8.5



DS-10



DS-11.5



DS-13

## 8 Pilot/Final Drill Stopper



PS-3.5

## 9 Fixture Driver

HD-IU-RL  
RD-IU-RS  
RD-IU-RL

## 10 Extension



WDE

## 11 1.2 Hex



HD-TW-L

(Bottom of the Case)  
Depth Gauge



DP

Torque Wrench




TW

# IU Stopper Drill Kit

# IU Stopper Drill Kit

Art No. : IU SD KIT



1 Point Drill



PD

2 Final Drill



Ø2.2  
TSD22070  
TSD22085  
TSD22100  
TSD22115  
TSD22130



F3.6  
ISD35070  
ISD30085  
ISD30100  
ISD30115  
ISD30130



F4.0  
ISD35070  
ISD35085  
ISD35100  
ISD35115  
ISD35130



F4.5  
ISD40070  
ISD40085  
ISD40100  
ISD40115  
ISD40130



F5.0  
ISD45070  
ISD45085  
ISD45100  
ISD45115  
ISD45130



F5.5  
ISD50070  
ISD50085  
ISD50100  
ISD50115  
ISD50130



F6.0  
ISD55070  
ISD55085  
ISD55100

3 F3.3

Lindemann Ø2.7 Twist Cortical



NLD-20



NSD-27



NPD-33

4 Cortical Drill



IUCD-36



IUCD-40




IUCD-45



IUCD-50



IUCD-55



IUCD-60

5 G/H Checker &amp; Path Finder



ICOGHCA



ICOGHC

6 Guide Pin



GPC

7 Extension



WDE

8 1.2 Hex



HD-TW-L

9 Fixture Driver

HD-IU-ML  
RD-IU-MLHD-IU-RL  
RD-IU-RS  
RD-IU-RL

10 Torque Wrench



TW

(Bottom of the Case)

Depth Gauge




DP



# Point Drill


- Initial drill for marking the implantation point easily on the cortical bone without slipping when making an implantation hole.

Art No.	Diameter(mm)	Edge Length(mm)
PD	1.5	9.5



# Stopper Final Drill

- Designed in accordance with the shape of the fixture to be implanted
- Easy identification of drills with color coding
- Adjusts drilling depth with a 1mm margin up to the fixed stopper part, including a laser-marked line




Fixture	F3.6	F4.0	F4.5	F5.0	F5.5	F6.0	Drill Length
7.0	ISD30070	ISD35070	ISD40070	ISD45070	ISD50070	ISD55070	30
8.5	ISD30085	ISD35085	ISD40085	ISD45085	ISD50085	ISD55085	31
10.0	ISD30100	ISD35100	ISD40100	ISD45100	ISD50100	ISD55100	32.5
11.5	ISD30115	ISD35115	ISD40115	ISD45115	ISD50115		32.5
13.0	ISD30130	ISD35130	ISD40130	ISD45130	ISD50130		34
Diameter(mm)	2.95	3.6	4.2	4.8	5.3	5.8	

Twist Stopper Drill

- Ø2.2 drill used for initial drilling to the length of the fixture to be implanted.
- Adjusts drilling depth with a 1 mm margin up to the fixed stopper part, including a laser-marked line

Art No.	Length(mm)
TSD22070	30
TSD22085	31
TSD22100	32.5
TSD22115	32.5
TSD22130	34



Stopper Cortical Drill

- Drill used for cortical bone removal in hard bone.
- Easy identification of drills with color coding

Art No.	Length(mm)	Diameter
IUCD-36	27.3	3.65
IUCD-40	27.3	4.15
IUCD-45	27.3	4.65
IUCD-50	27.3	5.15
IUCD-55	25.8	5.65
IUCD-60	25.8	6.15




## Narrow Lindemann Drill

\* Optional

- Lindemann drill designated for IU F3.3.
- Make a drill hole to determine the implantation depth and axis.
- Drilling as deep as the length of the fixture is required.
- Angle correction is available while drilling.

Art.No	Application
NLD-20	IU F3.3




## Narrow Straight Drill

\* Optional

- Straight drill designated for IU F3.3.
- Drilling as deep as the length of the fixture is required.
- The color coding on the upper part of drill is indicating the diameter and the main fixture that is being used.

Art.No	Application
NSD-27	IU F3.3




## Narrow Final Drill

\* Optional

- Final drill designated for IU F3.3.
- Drill for expanding the cortical bone.
- Use it after making the final holes on hard bones.


Art.No	Application
NPD-33	IU F3.3



Guide Pin

- Used to identify the direction and placement of the bone preparation after using the Lindemann Drill.

Art No.
GPC



Drill Extension

- An instrument to extend the length of the drill and the driver of the machine.
- Incomplete assembly may cause bending or fracture of the extension.
- The length can be extended by 15.4mm when the drill extension is assembled.


Art No.	Extendable Length(mm)
WDE	15.4



Gingiva Height Checker & Path Finder

- Check the height of the gingiva with the indicated color line.
- In the multiple case, assemble with the fixture that is placed firstly and confirm the precise implantation path.


Art No.	Application
IC0GHCA	MINI
IC0GHC	STANDARD



# Fixture Driver Machine Type


- Used when placing the fixture with a machine.
- Connections are divided into F3.3(Mini)/F4.0-6.0(Standard).
- Short and long drivers are provided for the convenience of the user.

Art No.	Length(mm)	Application	
HD-IU-MS	25	Mini	* Optional
HD-IU-ML	29		
HD-IU-RS	25	Standard	* Optional
HD-IU-RL	29		



# Fixture Driver Torque Wrench Type


- Used when placing the fixture with a torque wrench.
- Must check that it is completely assembled before use.  
Incomplete assembly may cause fracture due to heavy stress.
- Be careful, as excessive torque may cause damage to the internal hex.
- Short and long drivers are provided for the convenience of the user.



Art No.	Length(mm)	Application	
RD-IU-MS	13.2	Mini	* Optional
RD-IU-ML	18.8		
RD-IU-RS	13.2	Standard	* Optional
RD-IU-RL	18.8		

1.2 Hex Driver Torque Wrench type

- Can be used for both torque wrench & hand driver.
- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 45 Ncm. Exceeding the torque range may damage the tip.
- Various lengths are available for the convenience of the user.





Art.No	Length(mm)	
HD-TW-S	8mm	* Optional
HD-TW-L	13mm	
HD-TW-LL	18mm	* Optional

1.2 Hex Driver Machine Type

\* Optional

- Driver for Machine.
- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 45Ncm. Exceeding the torque range may damage the tip.
- Various lengths are available for the convenience of the user.


Art.No	Length(mm)
HD-LT-S	20
HD-LT-L	24
HD-LT-LL	28



Depth Gauge

- Instrument for measuring the drilling depth. (7-13mm)


Art No.
DP



Torque Wrench

- Assemble it to the adapter or driver to check torque during implantation or screw tightening.
- Apply torque by pulling the bar to the indicated torque value line the user wants to apply.
- Rotate the pivot hand of the torque wrench to change the torque direction.

Art No.
TW




KAR	ESP	WIRE	WISE II	WISE	PM Kit	WPK	WAGA	UT/IU Compact	IUT	IU Stopper Drill	IU Compact	IU Standard	IU Full
-----	-----	------	---------	------	--------	-----	------	---------------	-----	------------------	------------	-------------	---------





# Combined Surgical Kit for IU/UT Implant

Art No. : IUT KIT



1 Point Drill 2 Lindemann Drill



3 Bone Collecting Stopper



4 Bone Collecting Drill




5 Final Drill



6 Tap Drill




7 Fixture Driver



8 Guide Pin



9 G/H Checker & Path Finder



10 Extension




11 1.2 Hex



12 Torque Wrench



(Bottom of the Case)  
Depth Gauge



Tap Drill

- A dedicated tap for UT Implants.
- It is used to create a micro threaded channel in hard bone.
- Use it with the torque wrench after using the machine or the shank adapter.
- Tapping up to the bottom of the laser-marked line is recommended.
- Use 7mm(hard cortical bone), 10mm(thick and hard cortical bone)
- F3.6: Yellow, F4.3: Blue, F4.8: Pink, F5.3: Purple, F6.3: Brown



Fixture	F3.6	F4.3	F4.8	F5.3	F6.3
Art No.	TD36070N	TD43070N	TD48070N	TD53070N	TD63070N

KAR	ESP	WIRE	WISE II	WISE	PM Kit	WPK	WAGA	UT/IU Compact	IUT	IU Stopper Drill	IU Compact	IU Standard	IU Full
-----	-----	------	---------	------	--------	-----	------	---------------	-----	------------------	------------	-------------	---------


# UT Compact Kit

Compact Kit for UT Implant



# Compact Kit for UT Implant(F4.3-F4.8)


Art No. : UT BASIC KIT



- 1 Initial Drill      2 Pilot Drill      3 Step Drill      Tap Drill      Stopper
- LD-20-L      PD-30      SD43085      SD48085      TD43070N      TD48070N      S-43      S-48  
SD43100      SD48100      SD43115      SD48115
- 4 Initial Drill Stopper      5 Guide Pin      6 G/H Checker & Path Finder      7 Fixture Driver      8 Drill Extension
- DS-8.5      DS-10      DS-11.5      GPS      IC0GHCC      HD-IU-RL      RD-IU-RS      RD-IU-RL      WDE
- 9 1.2 Hex      10 Shank Adapter
- HD-TW-L      RA-HS-L      Torque Wrench      Depth Gauge      TW      DP
- (Bottom of the Case)

## Step Drill

- IT Implant drill by diameter and length.
- Color coding indicates the fixture's diameter.
- F3.6: Yellow, F4.3: Blue, F4.8: Red, F5.3: Purple, F6.3: Brown



Art. No.					Length(mm)
F3.8 <span style="color: orange;">●</span>	F4.3 <span style="color: blue;">●</span>	F4.8 <span style="color: red;">●</span>	F5.3 <span style="color: purple;">●</span>	F6.3 <span style="color: brown;">●</span>	
SD38070 *	SD43070 *	SD48070 *	SD53070 *	SD63070 *	7
SD38085 *	SD43085	SD48085	SD53085 *	SD63085 *	8.5
SD38100 *	SD43100	SD48100	SD53100 *	SD63100 *	10
SD38115 *	SD43115	SD48115	SD53115 *		11.5
SD38130 *	SD43130 *	SD48130 *	SD53130 *		13

\* : Optional

## Step Drill Stopper

- An instrument that helps the precise drilling to the depth of the fixture when using the step drill.
- By connecting the stopper, the length of step drill will be 0.5 mm longer, and additionally, by adding the 0.5 mm of Y-dimension, it will be 1 mm longer. Therefore, users must be careful to apply it to the lower anatomical structure.




Art. No.	Application	Drill Stop
S-38 *	F3.6	Step Drill 7.0 ~13.0mm (+1.0mm)
S-43	F4.3	
S-48	F4.8	
S-53 *	F5.3	

\* : Optional

Shank Adapter

- A connector that allows a machine driver to be assembled to a torque wrench.

Art. No.	Length	Type	
RA-HS-S	4.5	Short	* Optional
RA-HS-L	9.2	Long	





# WAGA Kit

Warantec Guide Kit




## WAGA KIT

- » **Warantec's Guide Kit can be used for the IU Implant System.**
  - Basic Composition : for the IU implant system
- » **A separate kit is required for fixtures with a wide diameter.**
  - Warantec's Guide Kit Plus : F5.5, F6.0
- » **Open type metal sleeve is convenient to use on the posterior area with a limit to opening the mouth.**
- » **Composed of 2 types of metal sleeves according to the diameter of guide hole.**
  - Regular Sleeve : F3.6, F4.0, F4.5, F5.0
  - Wide Sleeve : F5.5, F6.0
- » **Composed of Vertical/Lateral Anchor.**
  - The accuracy of the guide becomes excellent with the anchor, as the guide can be fixed stably.
- » **Sleevesless**
  - Composed of two library types (opened and closed) in accordance with the drill's diameter.

Applicable software : 3shape  exocad

## Sleeve

	Regular	Wide
Used Fixture	F3.6-F5.0	F5.5-F6.0
Inner Diameter	5.2	6.5
Outer Diameter	7.5	8.8
ART NO.	MSV-52	MSV-65




Warantec's guide system can be used with a sleeves or sleeveless options.

Based on your preference, you can choose either a sleeve or sleeveless design from our library.

# Warantec Guide Kit for IU Implant(F3.6-F5.0)

Art No. : WAGA KIT



1 Flattening Drill



2 Initial Drill



3 Step Drill



4 Tissue Punch



5 Bone Profile Drill & Screw



6 Final Drill



7 Fixture Driver



8 1.2 Hex Driver



11 Anchor Driver




10 Lateral Anchor Screw & Drill



9 Vertical Anchor for Standard(F4.0, F4.5, F5.0)




12 Torque Wrench



## Initial Drill

- Initial drill for marking the implantation point easily on the cortical bone without slipping when building an implantation hole.
- It can be used on offset by choosing the length between short and long.
- Recommended RPM : 800-1,000


Art No.	Diameter
GLD28	2.8



## Tissue Punch

- Used for gingival removal during flapless surgery.
- Recommended RPM : 800-1,000


Art No.	Diameter
GTP	4.5



## Flattening Drill


- Used to planarize the irregular bone area when the top of the alveolar bone is not flat, or it is narrow.
- Multiple cutting edge help to remove the bone stably without bouncing.
- Recommended RPM : 800-1,000

Art No.	Offset	Length	
GBF090	9	23.8	
GBF105	10.5	25.3	
GBF120	12	28.3	* Optional
GBF	No offset	27.5	* Optional



## Step Drill

- The saline flow design is applied to the part where the drill is guided to prevent bone overheating.
- Various lengths of drills are available for the convenience of the user.
- Recommended RPM : 800-1,000



Fixture	F3.6	F4.0	F4.5	F5.0	F5.5	F6.0	
Diameter	2.9	3.6	4.1	4.7	5.3	5.8	Length
Art No.	GSD30085	GSD35085	GSD40085	GSD45085	GSD50085	GSD55085	8.5
	GSD30100	GSD35100	GSD40100	GSD45100	GSD50100	GSD55100	10
	GSD30115	GSD35115	GSD40115	GSD45115	GSD50115	GSD55115	11.5
	GSD30130	GSD35130	GSD40130	GSD45130	GSD50130	GSD55130	13
	GSD30145	GSD35145	GSD40145	GSD45145	GSD50145		14.5
Sleeve	MSV-52				MSV-65		
Kit	WAGA Kit				WAGA Kit Plus		

## Final Drill

- Drill for expanding the cortical bone after using the bone collecting drill.
- Use it after making the final holes on hard bones.
- Use in consideration of the diameter of the fixture.
- Recommended RPM : 800-1,000




	F3.6	F4.0	F4.5	F5.0	F5.5	F6.0
Diameter	3.6	4.1	4.6	5.1	5.6	6.1
Art No.	GPD-35	GPD-40	GPD-45	GPD-50	GPD-55	GPD-60
Sleeve	MSV-52				MSV-65	
Kit	WAGA Kit				WAGA Kit Plus	

# Fixture Driver Machine Type


- Used when placing the fixture with a machine.
- Connections are divided into F3.3-F3.6(Mini)/F4.0-6.0(Standard).
- The outer hex indicator is applied to the driver to confirm the direction of the hex inside the fixture.

Art No.	Application
GFD-HM	Mini
GFD-HS	Standard



# Fixture Driver Torque Wrench Type


- Use it in case of placing the fixture by using a torque wrench.
- Must check if it is completely assembled during the usage. If it is assembled incompletely, it may cause fracture under the heavy stress.
- Be careful because the excessive torque may cause the damage of internal hex.
- The outer hex indicator is applied to the driver confirm the direction of the hex inside the fixture.



Art No.	Offset	Application	
GFD-RM120	12	Mini	
GFD-RS105	10.5	Standard	
GFD-RS090	9	Standard	* Optional

Bone Profile Drill


- A drill to remove bone after the placement of the fixture to prevent interference between the bone and the prosthesis.
- Assemble the guide screw to the fixture, and it allows the peripheral bone to be easily removed without damaging the fixture.
- Choose the drill to use by considering the type of the fixture and the size of the abutment.
- Recommended RPM : 800-1,000



Art No.	Diameter(mm)
BP-45	4.5
BP-55	5.5
BP-75	7.5

Guide Screw for Bone Profile Drill

- A guide screw that is assembled with the fixture and stops the bone profile drill vertically.
- Choose the screw to use by considering the type of the fixture.
- Assemble it by using the 1.2 hex driver with your hand.



Mini


Standard

Art No.	Application
IU-PGS-M	Mini
IU-PGS-S	Standard

## 1.2 Hex Driver Torque Wrench type


- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 35 - 45 Ncm.
- Exceeding the torque range may damage the tip.

Art No.	Hex Size
HD-TW-L	1.2



## Vertical Anchor

- A tool that fixes the guide by assembling vertically to the fixture.
- Usually used for edentulous jaw cases, and prevents vertical movement.
- Assemble it by using the vertical anchor driver with hand.



Art No.	Offset(mm)	Application	
GVA120M	12	F3.6	* Optional
GVA090S	9	F4.0-F5.0	
GVA105S	10.5	F4.0-F5.0	
GVA090W	9	F5.5-F6.0	* Optional
GVA105W	10.5	F5.5-F6.0	* Optional
GVA120W	12	F5.5-F6.0	* Optional



Lateral Anchor Drill

- A drill to create a hole on the lateral to place the Anchor Screw.
- Recommended RPM : 800-1,000


Art No.
GLAD



Lateral Anchor Screw & Pin

- The anchor screw that fixes the guide on the lateral.
- Assemble it by using the Lateral Anchor driver with your hand.


Art No.	Type	
GLAS	Screw	
GLAP	Pin	* Optional



Driver for Anchor


- A tool for assembling the Vertical Anchor & the Lateral Anchor Screw and Pin into the hole formed on the lateral.
- Assemble it by hand or using a torque wrench.

Art No.	Type
GAD-H	Machine
GAD-R	Torque Wrench



Torque Wrench

- Assemble it with the adapter or driver to check the torque during the placement of the fixture or tightening the screw.
- Shows torque by pulling the bar to the indicated torque value line that user wants to apply.
- Rotate the pivot hand of the torque wrench to change the torque direction.



Art No.

TW

Guide Bush

\* WAGA Kit Plus

- The axis guide tool for the case of using the small diameter drill on the wide sleeve.
- It is included in the Warantec's Guide Kit Plus.



Art No.


Inner Diameter

Outer Diameter

MGB-52

5.2

6.5




Guide Reamer Drill

\* Optional

- The reamer for polishing the hole after printing out the WAGA guide template with sleeveless.
- Recommended RPM for machine : 50-80 RPM
- Can be used with a Shank Adapter with hand.


Art No.	Diameter(mm)
GHR5	5.2
GHRW	6.5



Shank Adapter

\* Optional

- Required when using Reamer Drill.
- A connector that allows a driver for a handpiece to be connected to a torque wrench.



Short      Long

Art.No.	Length(mm)	Type
RA-HS-S	4.5	Short
RA-HS-L	9.2	Long

# WPK


Warantec Prosthetic Kit



IU Full
IU Standard
IU Compact
IU Stopper Drill
IUT
UT/IU Compact
WAGA
WPK
PM Kit
WISE
WISE II
WIRE
ESP
KAR

# Warantec Prosthetic Kit


Art No. : WPK KIT




1 MUA Adapter : Ø4.0




2 MUA Adapter : Ø5.0




3 SB Adapter for only IT System




4 TOP Adapter for only IT System



5 1.2 Hex Driver




6 Shank Adapter



7 Extension Smart Driver Extension




8 G/H Checker & Path Finder



## Hex Adapter Machine Type

- Adapter for machine.
- Used to assemble the multi-unit abutment to the fixture.


Art No.	Diameter(mm)	Length(mm)	Application
HA-HHA	4.0	14	MUA 4.0
HA-HH	5.0	18	SUA 5.0 & MUA .0



## Hex Adapter Torque Wrench Type

- Can be used for both torque wrench & hand adapter.
- Used to assemble the multi-unit abutment to the fixture.
- Various lengths of adapters are available for the convenience of the user.


Art No.	Diameter(mm)	Length(mm)	Application
RA-HHA-S	4.0	6.6	SUA 4.0
RA-HHA-L	4.0	11.1	
RA-HHA-LL	4.0	16	
RA-HH-S	5.0	6.6	SUA 5.0 & MUA 5.0
RA-HH-L	5.0	11.1	
RA-HH-LL	5.0	16	



## 1.2 Hex Driver Machine Type

- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 35 - 45 Ncm.
- Exceeding the torque range may damage the tip.
- Various lengths of drivers are available for the convenience of the user.


Art.No	Length(mm)
HD-LT-S	20
HD-LT-L	24
HD-LT-LL	28



## 1.2 Hex Driver Torque Wrench type

- Used for cover screw, abutment screw, and healing abutment.
- The hex tip is designed to tolerate the maximum torque of 35 - 45 Ncm.
- Exceeding the torque range may damage the tip.
- Various lengths of drivers are available for the convenience of the user.


Art No.	Length(mm)	
HD-TW-S	8	
HD-TW-L	13	
HD-TW-LL	18	* Optional



## Shank Adapter

- Required when using Reamer Drill.
- A connector that allows a driver for a handpiece to be connected to a torque wrench.

Art No.	Length(mm)
RA-HS-S	4.5
RA-HS-L	9.2



Drill Extension

- An instrument to extend the length of the drill and the driver of the machine.
- Incomplete assembly may cause bending or fracture of the extension.
- The length can be extended by 15.4 mm when the drill extension is assembled.


Art No.	Extendable Length (mm)
WDE	15.4



Gingiva Height Checker & Path Finder

- Check the height of the gingiva with the indicated color line.
- In the multiple case, assemble with the fixture that is placed firstly and confirm the precise implantation path.



Art No.	Application
IC0GHCA	Mini
IC0GHC	Standard





**SDE(Smart Driver Extension)**

- 1.2 hex driver extension.
- Used when access is difficult due to the long Anatomical Crown Length.



Art No.	Anatomical Crown Length(mm)
SDE-S	20
SDE-L	23
SDE-LL	26


# PM Kit

## Prosthetics Multi Kit



# Prosthetics Multi Kit

Art No. : PM KIT



- ① Adapter for Ø4.0 Multi-unit Abutment(Torque Wrench)      Driver for Ø5.0 Multi-unit Abutment(Torque Wrench)



- ② KAR Lead Driver      KAR Hook Housing




- ③ ESP Shaft      ESP Holder



- ④ Elastic Gasket




- ⑤ G/H Checker & Path Finder      ⑥ Smart Extension      ⑦ 1.2 Hex



## Hex Adapter Torque Wrench Type

- Can be used for use both torque wrench & hand adapter.
- Used to assemble the multi-unit abutment to the fixture.
- Various lengths of adapters are available for the convenience of the user.






Art No.	Diameter(mm)	Length(mm)	Application	
RA-HHA-S	4.0	6.6	MUA 4.0	* Optional
RA-HHA-L	4.0	11.1		
RA-HHA-LL	4.0	16		* Optional

Art No.	Diameter(mm)	Length(mm)	Application	
RA-HH-S	5.0	6.6	SUA 5.0 & MUA 5.0	* Optional
RA-HH-L	5.0	11.1		
RA-HH-LL	5.0	16		* Optional

## SDE(Smart Driver Extension)

- 1.2 hex extension driver.
- Used when access is difficult due to the long Anatomical Crown Length.




SDE-S SDE-L SDE-LL

Art No.	Anatomical Crown Length(mm)	
SDE-S	20	* Optional
SDE-L	23	
SDE-LL	26	* Optional

## 1.2 Hex Driver Torque Wrench type

- 1.2 hex driver.
- Assemble only by hand force using the ESP shaft.
- Short, long, and extra long drivers are provided for the convenience of the operator.




Art.No	Length(mm)	
HD-TW-S	8	* Optional
HD-TW-L	13	
HD-TW-LL	18	* Optional

## Gingiva Height Checker & Path Finder


- Check the height of the gingiva with the indicated color line.
- In the multiple case, assemble with the fixture that is placed firstly and confirm the precise implantation path.

Art No.	Application
IC0GHCA	Mini
IC0GHC	Standard



## ESP Shaft

- Easy handling as the abutment and the crown can be held simultaneously.
- Head part can be assembled with a 1.2 hex driver.
- Added the 1.2 hex design on the bottom so the screw can be assembled using a shaft.
- Applied the safety design so that the head part gets fractured if the torque value exceeds the appropriate (25N) value.
- Recommended tightening torque : 15Ncm



Art No.	Length(mm)
ESP-S1825	25

## ESP Holder

- Fix the abutment and crown by assembling them with the shaft.
- Use it for Try-in & out.



ESP-H

## Elastic Gasket

- Hold the abutment by assembling with the shaft.
- Applied the design that can easily remove the remaining cement inside the crown.
- It is a disposable product. Reusing the product may cause a technical problem.
- Easy replacement with gasket holder.



ESP-EG18




Gasket Holder (ORH)

Please refer to the ESP category for the instructions and optional products.

## Lead Driver


- It can be used to remove the abutment when the assembled part is fractured with hook housing.
- A tool that is inserted by penetrating the Hook Housing.
- Blocks contraction to maintain the grip force of Hook Housing that holds the abutment, and at the same time, disassembles the abutment from the fixture by pushing the bottom of the fixture with the turning force that is generated from engaging it to the Hook Housing.
- Use a torque wrench.
- Recommended : Less than 50 Ncm



Art No.	Length(mm)
KAR-LD20425	41.2

## Hook Housing


- A tool that fixes vertically on the tip of hex or non-hex of the abutment that is engaged with the fixture by penetrating the screw hole of the abutment.
- It can be used to remove the abutment when the assembled part is fractured with lead driver.



Art No.	Length(mm)
KAR-HH20425	27.8

## Housing Holder

- A tool that holds the Hook Housing so that it does not rotate when engaging the lead driver into hook housing.
- The round tip on the other side can be used when the operator uses a destructive method to remove the abutment.



Art No.
KAR-HH

Please refer to the KAR category for the instructions and optional products.

# WISE Kit

Warantec Implant's Sinus Elevation Kit



IU Full
IU Standard
IU Compact
IU Stopper Drill
IUT
UT/IU Compact
WAGA
WPK
PM Kit
WISE
WISE II
WIRE
ESP
KAR




# Sinus Combination Kit for Crestal & Lateral Approach

Art No. : WISE KIT



WISE KIT CASE



15 Aqua Silicone Tube



14 Bonewell



16 Flexible Sinus Curette



7 Lateral Stopper



8 Sinus Reamer



9 Core Drill



10 Sinus Diamond Reamer



11 Diamond Core Reamer



12 Side Drill



13 Depth Gauge



## Crestal Approach Components

Sinus Point Drill



Ø1.8    Ø2.3

Crestal Reamer



Ø2.8    Ø3.3    Ø3.8    Ø4.2

Crestal Stopper



Crestal Diamond Reamer



Ø2.8

Aqua Lift Tap



Ø2.8/3.8    Ø3.3/4.2

Aqua Hand Adapter



Depth Gauge




Aqua Silicone Tube



Bonewell



### How to Use



- Select the safe and appropriate crestal stopper according to the bone thickness on the patient's x-ray.

- The Reamer Drill must be used with the selected drill stopper for accurate and safe drilling.

- Use a Sinus Reamer with an appropriate diameter (2.8, 3.3, 3.8, 4.2) according to the size of the implant and bone density. (800-1,000RPM)

- Use diamond bur to confirm the final detachment of the sinus membrane and to shape any irregularity of the sinus floor.

- Lift the sinus membrane safely with the water pressure using the Aqua Lift Tap.

## Lateral Approach Components

Lateral Reamer



Lateral Core Drill




Diamond Reamer &amp; Core Reamer




Lateral Side Drill



Flexible Sinus Curette




Lateral Stopper




### How to Use

①




②



- Select the safe and appropriate lateral stopper according to the bone thickness on the patient's x-ray
- The Core Drill or Diamond Core Reamer must be used with the selected drill stopper for accurate and safe drilling.
- Use the Core Drill (Ø7.0) or Diamond Core Reamer (Ø7.0) depending on the thickness and density of the bone for safe drilling. (800-1,000RPM)
- The sinus membrane can be safely detached by using the Easy Sinus Flexible Curette.


③




④



①



②



- Select the safe and appropriate lateral stopper according to the bone thickness on the patient's x-ray.
- The Sinus Reamer or Diamond Reamer must be used with the selected drill stopper for accurate and safe drilling.
- For safe drilling (800-1,000RPM), choose Sinus Reamer (Ø8.0, Ø6.5/8.0) or Diamond Reamer (Ø6.5/8.0) according to the thickness and density of the bone.
- The sinus membrane can be safely detached by using the Easy Sinus Flexible Curette.


③




④



⑤



⑥



# WISE Kit II



Warantec Implant's Sinus Elevation Kit II



# Warantec Implant's Sinus Elevation Kit II

Art No. : WISE Kit II

WISE KIT II is a sinus master kit that can proceed Crestal approach & Lateral approach with the pop-up drill and Aqua lift system configurations that make sinus elevation the easiest, simplest, and safest.







## ① Pop-up Initial Drill

- Drill forming a hole on the implantation point.
- Forming a Guide hole for the pop-up bar of pop-up drill.
- Secure 2mm of residual bone for the safe elevation.
- Recommended RPM : 500 - 800


Art No.
OSD-32T



## ② Pop-up Drill

- Safely elevate the membrane with the pop-up function.
- Drill it after inserting the pop-up bar into the hole that formed by initial drill.
- Engage the stopper for the safe drilling.
- Recommended RPM: 200 - 500

Art No.	Diameter (Ø)
AD-32	Ø3.2
AD-38	Ø3.8



## ③ Pop-up Drill Stopper

- Drill stopper that can safely control the straight moving of the drill with the spring's elasticity.
- The distance from the blade at the end to the stopper is 3mm when it is engaged with the pop-up drill.
- While drilling, in case of the spring is completely compressed, it can stop 5mm.


Art No.
AD-ST



④ Pop-up Drill Remover

- Tool for disassemble the drill and the shank parts.
- After using the drill, make sure to disassemble and clean the drill, then reassemble and sterilize it.

Art No.
ADR



⑤ Bone Pusher


- This tool is utilized for grafting collected bone or bone material after elevating the sinus membrane.
- It functions as an osteotome when elevating the sinus membrane using a mallet.

Art No.
BP-35C




⑥ Aqua Lift System

- Secure the aqua lift tab in place by utilizing the aqua hand adapter within the hole formed on the maxillary sinus bone.
- Connect the aqua lift tab to the syringe tube and inject saline into the hole to elevate the maxillary sinus.
- Insert the syringe spring into the syringe to facilitate convenient water pressure adjustment.




Art.No.	Product
ALT-28	Aqua Lift Tap
AQUA-HA	Aqua Hand Adapter
AST-30	Syringe Tube
S-Spring	Syringe Spring



# Lateral Side Drill

- The Drill for extending the hole created by the Pop-up Drill.
- Recommended RPM: 800 - 1,000



Art No.
SID-30



# Sinus Curette

- Utilized for safely elevating the maxillary sinus membrane or performing bone graft.

Art No.
7091
7094



# Depth Guage

- Use the Depth Gauge to inspect for any perforations and determine the drilling depth of the maxillary cortical bone.
- Additionally, it is used for elevating bone grafts following maxillary sinus elevation.

Art No.	
DG-01	* Optional



# Bonewell



Art No.	
AT-02	* Optional


# Aqua silicone tube Ø3.0



Art No.	
5702	* Optional

# Sinus Curette

Art No.	
7092	* Optional
7093	* Optional






## SURGICAL PROCEDURE FOR CRESTAL APPROACH


### 1. Marking (Use Pop-up Initial Drill)

- Utilize the pop-up initial drill to mark the precise point for implant placement.
- For the secure elevation of the pop-up drill, drill until at least 2mm of residual bone remains.



### 2. Bone Lifting (Use Pop-up Drill)

- Insert the pop-up bar of the pop-up drill into the designated hole and proceed towards the sinus while drilling.
- As the drill approaches the sinus, the pop-up bar's vertical pressure lifts the residual bone and elevates the membrane.
- To avoid potential damage to the sinus membrane while drilling, we recommend to use the spring stopper. (Recommended RPM: 200 - 500)



※Tip : If the residual bone in contact with the sinus is hard, making the final bone lifting challenging, you can utilize the Osteotome technique with the assistance of the Bone Pusher to lift the residual bone.


### 3. Sinus Membrane Lifting (Use Aqua Lift System)

- Following the bone lifting, utilize the adapter to connect the aqua tip with the syringe.
- Elevate the sinus membrane using water pressure while injecting saline through the syringe.
- For a smoother elevation process and to avoid sudden water pressure, insert a syringe spring into the syringe.



### 4. Implant Placement(with GBR or without GBR)


- Fill the elevated space with bone material for the subsequent implantation, or proceed the treatment without the bone graft.



## SURGICAL PROCEDURE FOR LATERAL APPROACH


### 1. **Marking** (Use Pop-up Initial Drill)

- Utilize the pop-up initial drill to mark the precise point for implant placement.
- For the secure elevation of the pop-up drill, drill until at least 2mm of residual bone remains.



### 2. **Bone Lifting & Cutting** (Use Pop-up Drill and Side Drill)

- Insert the pop-up bar of the pop-up drill into the designated hole and proceed towards the sinus while drilling.
- As the drill approaches the sinus, the pop-up bar's vertical pressure lifts the residual bone and elevates the membrane.
- To avoid potential damage to the sinus membrane while drilling, we recommend to use the spring stopper. (Recommended RPM: 200 - 500)
- Extend the hole for elevation using the side drill.




### 3. **Sinus Membrane Lifting** (Use Sinus Curette)

- After elevating the bone, carefully lift the membrane using a sinus curette.



### 4. **GBR & implantation**

- Fill the elevated space with the bone material.
- After completing the bone graft, proceed with the implantation.




# WIRE Kit

Warantec Implant's Rescue for Emergency Kit



# Warantec Implant's Rescue for Emergency

Art No. : WIRE KIT



1 Fixture Removal Tool



FR-IU

2 Trephine Drill



TPD-35 TPD-40 TPD-45 TPD-55

3 Abutment Removal Drill & Tool



SD-AR-16



RD-AR-16

4 Dual Abutment Removal Tool



RD-AR-S



RD-AR-L

5 Screw Removal Tool



HD-SR-CM

6 Screw Removal Drill & Tool



SD-SR-08

HD-SR-08

SD-SR-12

HD-SR-12

7 IU Screw Removal Drill Guide



UG-SR-HA



UG-SR-H

8 Removal Guide Holder



RGH

9 Remover Torque Wrench



RTW


The **WIRE KIT** is a comprehensive removal kit that enables removing implants and debris remaining on implants in case of failure of prosthetics such as a screw or abutment fracture.

**WIRE**

## Fixture Remover

- Implant fixture removal tool.
- When it is necessary to remove the implant due to peri-implantitis or other causes.
- Only the implant can be removed without damaging the surrounding bones.
- Easy and convenient removal of the internal fixture by reverse rotation.

❖ It is a disposable product. Reusing the product may cause a problem in its function.



FR-IU



1  
Tighten the fixture remover to the fixture to be removed by rotating it in the reverse direction.




2  
Maintain the fixture remover vertically, and use the remover torque wrench (RTW) to remove the fixture by rotating the fixture remover in the reverse direction.

## Trephine Drill

- Drill for fixture removal.
- Check whether the size fits the fixture to be removed before use.
- Form a path on the bone surface by the reverse rotation, and then drill it with normal rotation.
- Recommended RPM : 800-1,200 RPM

Art No.	A	B
TPD-35	3.5	4.1
TPD-40	4.0	4.6
TPD-45	4.5	5.3
TPD-55	5.5	6.3



1 Assemble the trephine drill that corresponds to the fixture's diameter to be removed to the handpiece. Form a path on the bone surface by the reverse rotation.



2 Perform the normal rotation drilling up to the implanted fixture depth to separate the fixture to be removed from the surrounding bone.


\* Increase the amount of irrigation sufficiently to reduce bone heating.



3 In case of broken fixture, tighten the fixture remover by reverse rotation and remove the fixture.

## Abutment Removal Drill

- Removal drill for neck fractured solid abutment.
- Drill for forming a hole on the cross section of the fractured solid abutment with reverse rotation.
- Recommended RPM : 1,500 RPM
- Package Component : SD-AR-16



SD-AR-16

## Solid Abutment Removal Tool


- Removal tool for the neck fractured solid abutment.
- Remove by tightening the removal tool by reverse rotation on the hole formed on the cross section of the fractured solid abutment.
- Package Component : RD-AR-16



RD-AR-16



1  
Remove the soft tissue to the extent that the Abutment Removal Drill can be used.




2  
Proceed with the reverse rotation drilling on the cross section of the fractured solid abutment with the SD-AR-16 drill. Proceed with the drilling at the speed of 1,500 RPM or higher. Sufficient irrigation will be required as cutting heat will be generated by the drilling speed. (\*minimum cutting depth of 2mm)

☆TIP

If the remaining piece of the solid abutment comes out while proceeding with the reverse rotation drilling with the SD-AR-16 drill, remove the remaining abutment with a tweezer.



3  
Insert the RD-AR-16 into the formed hole.



4  
Reverse rotate the RD-AR-16 to remove the fractured abutment.


## Dual Abutment Removal Tool

- Removal tool for the neck fractured dual abutment.
- When the dual abutment is fractured on the inside of the fixture.
- When the hex of the abutment screw is damaged.
- Remove it by reverse rotation from inside the fractured dual abutment or on the screw hex.

❖ It is a disposable product. Reusing the product may cause a problem in its function.



RD-AR-S




RD-AR-L


### Neck Fractured Dual abutment Removal Procedure



1  
Insert into the hole of the fractured abutment with the Dual Abutment Removal Tool.



2  
Reverse-rotate the Dual Abutment Removal Tool until it is fitted with the fractured abutment.




3  
If it is fitted adequately, use forceps or any other tools to shake it and remove it.

### Removal Procedure for the Damaged Abutment Screw Hex



1  
When the 1.2 hex driver cannot be assembled because the hex part of the abutment screw is damaged.



2  
Insert the Dual Abutment Removal Tool into the damaged hex.




3  
Reverse-rotate the Dual Abutment Removal Tool to remove the damaged abutment screw.




Screw Removal Tool


- Removal tool for the fractured abutment screw.
- In case the screw is fractured inside the fixture.
- Going through the incomplete cross section on the upper part of the fractured screw to remove it.
- Recommended RPM : Less than 50 RPM



Assemble the HD-SR-CM with a handpiece, and insert it into the fractured screw.



Reverse-rotate it with a speed of less than 50RPM to remove the screw.




HD-SR-CM

Art No.	Length (mm)
HD-SR-CM	24

Screw Removal Drill

- Removal drill for the fractured abutment screw.
- In case the screw is fractured inside the fixture.
- A drill to form a hole in the incomplete cross section on the upper part of the fractured screw.
- The operator can identify the depth of the deletion through the laser marking on the upper part.
- Recommended RPM : 1,500 RPM

Art No.	Diameter (mm)
SD-SR-08	0.8
SD-SR-12	1.2



Screw Removal Tool


- Removal tool for the fractured abutment screw.
- Reverse-rotate it on the hole formed on the cross section of the fractured screw to remove it.

Art No.	Diameter (mm)
HD-SR-08	0.8
HD-SR-12	1.2




# Screw Removal Drill Guide


- A guide for removing the abutment screw.
- Supports to form a hole in the middle of the cross section of the fractured screw.




Art No.	Hex Dimension	Length (mm)	Application
UG-SR-HA	2.1	12.4	Mini
UG-SR-H	2.5	12.4	Standard



Assemble the Guide with the Guide Holder, and insert it into the fixture.



Hold the handle of the holder, and insert the Screw Removal Drill into the Guide.



Reverse-rotate the Screw Removal Drill to form a hole on the fractured screw. Proceed with the drilling at the speed of 1,500RPM or higher. Sufficient irrigation is required as cutting heat will be generated by the drilling speed. (\*minimum cutting depth of 2mm)



Insert the Screw Removal Tool into the hole of the fractured screw.



Reverse-rotate the Screw Removal Tool to remove the fractured abutment screw. (\*Remove it with the speed of less than 50 RPM while applying moderate vertical pressure.)

# Screw Removal Guide Holder

- Holder for screw removal guide.
- A guide for safely using the Screw Removal Drill in the oral cavity.



RGH

Art No.

RGH

# Remover Torque Wrench

- Torque wrench for remover.
- 0-100Ncm : Gauge for removing the abutment and screw.
- 0-400Ncm : Gauge for removing the fixture.



RTW

Art No.

RTW

# ESP Kit

Easy Screwmentation Positioner



# Easy Screwmentation Positioner

Art No. : ESP KIT

① ESP Shaft(2ea)

Screw Hole  $\varnothing 2.35 - 2.65$



ESP-S1810



ESP-S1815



ESP-S1820



ESP-S1825



ESP-S1830

\*Optional

② ESP Shaft(2ea)

Screw Hole  $\varnothing 2.10 - 2.30$



ESP-S1610




ESP-S1615



ESP-S1620



ESP-S1625



③ 1.2 Hex



HD-TW-S

④ ESP Holder(6 pcs)



ESP-H




O-ring

⑤ Elastic Gasket(3 sets)



ESP-EG18



ESP-EG16

※ Optional





ESP can not only use the SCRP bridge crown intraorally more easily, quickly, and safely, but also enables humidity-free bonding & cementation outside of the mouth. ESP with these characteristics can be a perfect tool for the modelless prosthetic work through digital. It is a tool for the rapid & perfect clean-up of remaining cement. Also, these advantages enable the intraoral cementation, compensating for the inaccuracy generated at the stages of impression taking and the final prosthesis assembly.




## ESP Shaft

- Easy handling as the abutment and the crown can be held simultaneously.
- Head part can be used as a hex driver.
- Added the 1.2 hex design on the bottom so the screw can be assembled with a shaft.
- Applied the safety design that makes the head part gets fractured if the torque value exceeds the appropriate (25N) value.
- Assemble by hand using a 1.2 hex driver.
- Recommended tightening torque : 15Ncm




Art. No.	Package Components	Applicable Screw Hole	Length
ESP-S1810		Ø2.35 - 2.65	10
ESP-S1815			15
ESP-S1820			20
ESP-S1825			25
ESP-S1830			30

### \* Optional


Art. No.	Package Components	Applicable Screw Hole	Length
ESP-S1610		Ø2.10 - 2.30	10
ESP-S1615			15
ESP-S1620			20
ESP-S1625			25

## Elastic Gasket


- Fix the abutment by assembling with the shaft.
- Applied the design that can easily remove the remaining cement inside the crown.
- It is a disposable product. Reusing the product may cause a technical problem.
- Easy replacement with gasket holder.




### How to Assemble the Elastic Gasket to a Shaft




Hold the shaft and the gasket holder.




Engage the hex part of a shaft to the cylinder part of the gasket holder.



Move one elastic gasket from the gasket holder to the shaft.



Check if the elastic gasket is placed on the shaft correctly.



Detach the shaft from the gasket holder.

## ESP Holder & Hex Driver

### ESP Holder

- It fixes the abutment and crown by being combined with the a shaft.
- Use it for Try-in & out.




ESP-H



### 1.2 Hex Driver

- Use it only with hand force by assembling it with the ESP shaft.
- Lengths of short, long, and extra long drivers are provided for the convenience of the user.
- Package components : hex driver




Art No.	Diameter(mm)	Length(mm)
HD-TW-L	1.2	13




Bone Level Protocol


Try-in



Try-out




Cementation




## Tissue Level Protocol


### Try-in




### Try-out




### Cementation



ESP Try-in



Sit the crown on the abutment using a shaft, and grab the abutment and crown as one piece by using the holder.



The operator can easily, quickly, and safely try-in the combined prosthesis through ESP.

Disengage the ESP shaft and holder.

Engage the abutment screw and driver from the outside of the mouth, then assemble it to the intraoral prosthesis.


## ESP Try-out



Tighten the prosthetic screw with a driver, then rotate it counterclockwise to disassemble.




Engage the ESP shaft and holder from the outside of the mouth, then assemble it on the intraoral prosthesis.




The operator can easily, quickly, and safely try-out the combined prosthesis through ESP.

## ESP Cementation



Engage the ESP to the abutment and proceed with the primer & bonding from outside of the mouth in humidity-free condition.



Combine the cemented prosthesis in humidity-free condition.

Remove the remaining cement in the marginal zone, disassemble the ESP shaft from the cemented prosthesis vertically, and thoroughly remove the remaining cement inside the prosthesis hole.

## ESP Cementation



Assemble the cemented prosthesis to the placed fixture.

Engage the abutment screw and driver from the outside of the mouth, then assemble it to the intraoral prosthesis.



Harden the cement after biting.


# KAR Kit

Kinetic Abutment Remover



# Kinetic Abutment Remover

Art No. : KAR KIT



- ① KAR Lead Driver
  - KAR-LD20325
  - KAR-LD20375
  - KAR-LD20425
  - KAR-LD20475
- ② KAR Hook Housing
  - KAR-HH20325
  - KAR-HH20375
  - KAR-HH20425
  - KAR-HH20475
- ③ KAR Housing Holder
  - KAR-HH




KAR is a tool that helps disengaging the abutment from the internal conical joint fixture when it is stuck due to friction without damaging the prosthesis through the controlled axial load.

Remnants of the fractured abutment can also be easily removed by using KAR.

**KAR**

## Lead Driver

- It can be used to remove the abutment when the assembled part is fractured with hook housing.
- A tool that is inserted by penetrating the Hook Housing.
- Blocks contraction to maintain the grip force of Hook Housing that holds the abutment, and at the same time, disassembles the abutment from the fixture by pushing the bottom of the fixture with the turning force that is generated from engaging it to the Hook Housing.
- Recommended : Less than 50 RPM




Art No.	Package Components	Length	Applicable Screw Hole Size
KAR-LD20325		32.5	1.95~2.15
KAR-LD20375		37.5	
KAR-LD20425		42.5	
KAR-LD20475		47.5	

### \* Optional

Art No.	Package Components	Length	Applicable Screw Hole Size
KAR-LD16425		42.5	1.8~1.95
KAR-LD16475		47.5	

# Hook Housing

- An instrument that goes through the screw hole of abutments to vertically reach the top part of the hex or non-hex structure in the abutment assembled to the fixture.
- It can be used to remove the abutment when the assembled part is fractured.



Art No.	Package Components	Length	Applicable Screw Hole Size
KAR-HH20325		19.1	1.95~2.15
KAR-HH20375		24.1	
KAR-HH20425		29.1	
KAR-HH20475		34.1	

**\* Optional**

Art No.	Package Components	Length	Applicable Screw Hole Size
KAR-HH16425		29.1	1.8~1.95
KAR-HH16475		34.1	

## Housing Holder

- A tool that holds the Hook Housing so that it does not rotate when engaging the lead driver into hook housing.
- The round tip on the other side can be used when the operator uses a destructive method to remove the abutment.




Art No.

KAR-HH

## Torque Wrench

\* Optional

- Used to apply torque by engaging it to the Lead Driver.
- Apply the torque by pulling the bar up to the indicated torque value line set by user.
- Rotate the pivot hand of the torque wrench to change the torque direction.




Art No.

TW

# Try-in Abutment Remover

## KAR Length Guide




Anatomical Crown Length

Usage Range Screw Hole	Anatomical Crown Length <sup>mm</sup>			
	6	11	16	21
1.8 ~ 1.95 mm			KAR-LD16425	KAR-LD16475
1.8 ~ 1.95 mm			KAR-HH16425	KAR-HH16475
1.95 ~ 2.15 mm	KAR-LD20325	KAR-LD20375	KAR-LD20425	KAR-LD20475
1.95 ~ 2.15 mm	KAR-HH20325	KAR-HH20375	KAR-HH20425	KAR-HH20475

1.8 ~ 1.95 mm	Hole size for M1.8 screw products
1.95 ~ 2.15 mm	Hole size for M2.0 screw products

## Anatomical Crown Length Try-In



If the color coding is not visible, **it cannot be operated.**


If the color coding is visible, **it can be operated.**

The color coding is visible, but it may interfere the adjacent tooth, **so choose the length by considering the situation.**

## Try-in Abutment Remover




Hold the Hook Housing with a Housing Holder, and rotate the Lead Driver with your hand until there is resistance.



Engage the Torque Wrench to the Lead Driver, and rotate it to insert it into the Hook Housing. It will push the bottom of the screw hole of the fixture, disassembling the abutment with a friction-cracking sound.

Try-in Abutment Remover




Reverse-rotate the Lead Driver to disengage, gently hold the head of the Hook Housing with hand, and remove it vertically from the abutment.



KAR	ESP	WIRE	WISE KIT II	WISE	PM Kit	WPK	WAGA KIT	UT/IU Compact	IUT	IU Stopper Drill	IU Compact	IU Standard	IU Full
-----	-----	------	-------------	------	--------	-----	----------	---------------	-----	------------------	------------	-------------	---------

# Fractured Abutment Remover


## ■ Fractured Abutment Remover



## Fractured Abutment Remover



Engage the Torque Wrench to the Lead Driver, and rotate it to insert it into the Hook Housing. It will push the bottom of the screw hole of the fixture, disassembling the remaining abutment with a friction-cracking sound.



Reverse-rotate the Lead Driver to disengage, gently hold the head of the Hook Housing with your hand, and remove it vertically from the remaining abutment.



## Fractured Abutment Remover




### Special Tip

A destructive method when it cannot be removed because of the strong friction:

- Insert the straight bar of the Housing Holder into the abutment and shake it left and right to break the friction.

## Operation Principle



Even if the system does not have enough space for the Hook Housing to engage to the bottom part, the Hook Housing will grab the internal wall of the screw hole and remove the abutment.

# WARANTEC

## New Products

**WA Motor**



**WA Bone**



**Bone Dozer**



**WA Mem**



### WA Motor

- Reusable irrigation tube after the sterilization.
- Auto calibration function.
- Powerful force up to 70 Ncm(20 :1 gear handpiece)
- Function to stop automatically when exceed the set torque value.
- Large LCD display for various information.
- Show torque and RPM in real-time.
- Optic motor.



Art. No.	Product
WA Motor	Implant Engine

### Thread Cutting Function

Thread cutting is available in Tapping mode among 6 basic programs.

[The thread cutter function supports the surgery when makes tap for implants in to dense (D1) bone.]



### WA Bone


- Bovine bone. (Cancellous 100%)
- High stability implantation through the natural bone mineral materials.
- Uniform particle formation that consisted of cancellous bone.
- Porous structure which is similar to cancellous bone of human.
- Accelerate the formation and growth of the new bone
- High biocompatibility reduces the inflammation and immune response.
- Vial Type.
- Particle size : 0.2-1.0 mm



Double-sealed  
Safety package



Wide internal surface area



Stable Porous Structure



Minimizes vitrification



Surface with the high bone  
conduction

Art. No.	Product
01 - 0210	WA Bone 0.1g
25 - 0210	WA Bone 0.25g
05 - 0210	WA Bone 0.5g
10 - 0210	WA Bone 1.0g

### WA Bone thesis statement



WA Bone and PRP  
Techniques




Sinus Bone Graft



Crestal(Alveolar Bone Graft)


### Easy Bone Dozer

- A disposable scraper for collecting autogenous bone.
- It has a curved structure in front so that it is easier for lateral access intraorally.
- The blade rotates 360°, so if you rotate it 90° to use it, you can use 4 new blades.
- User can check the amount of autologous bone collected in real time because the product is transparent.



Art. No.	Product
3500	Easy Bone Dozer

※ Please sterilize the product before use because it is packed in a non-sterile package.



Rotate it with 1.2 hex driver.



If the blade becomes dull, rotate it by 90° and use it upto 4 times.



Start collecting bones.




Bone collection completed.

### Indication

- For augmentation of bone surrounding implants placed in extraction sockets
- Ridge augmentation / reconstruction
- Filling of bone defects after root resection, regained teeth, and etc.
- Guided bone regeneration in dehiscence defects
- Guided tissue regeneration in periodontal defects
- Made in Korea / Immediately available / CE and FDA to be done by 2Q of 2024




### Features



Dense collagen pore structure




Easy manipulation (Thin & Flexible)  
Size 15 X 30 X 0.3mm (One size)



Over 99% purity collagen









**Address** 411~412, 474, Dunchon-dearo, Jungwon-gu, Seongnma-si, Gyeonggi-do, Republic of Korea  
**Tel** +82-2-3675-5851 **Fax** +82-2-3675-5853 **E-mail** [info@warantec.com](mailto:info@warantec.com) **www.WARANTEC.com**